


ENDOTRACHEAL TUBE INTRODUCER (BOUGIE)



Dynarex Resp-O2[™] Endotracheal Tube Introducers (Bougies) are premium medical devices designed to effectively guide and support endotracheal tubes during difficult intubation procedures. The sterile Introducer is disposable and comes individually wrapped for optimal hygiene.

- Designed to guide and support endotracheal tubes into patients' airways
- Ideal for difficult intubations and "blind" intubations
- Crafted from high-quality plastic materials for optimal durability
- Features clearly visible depth markings for improved safety
- Available in Child and Adult sizes
- Disposable and individually wrapped for optimal hygiene
- Not made with natural rubber latex



10 Fr15 FrBougieBougieREORDER NO. 36260REORDER NO. 36261PEDIATRICADULTPACKAGE QUANTITY
10/BoxPACKAGE QUANTITY
10/Box

Contact your local sales representative or call to find your nearest distributor.



L 1.888 DYNAREX

WE CARE LIKE FAMILY

NAREX 💮 DYNAREX.COM

• 10 GLENSHAW STREET | ORANGEBURG, NY 10962

SEBURG, NY 10962 ®2022 Dynarex Corporation

R220623