




LARYNGEAL MASK AIRWAYS

Silicone


Dynarex Resp-O₂™ Laryngeal Mask Airways (LMA) is a premium-quality respiratory device designed to help open and maintain airways in unconscious patients. The Laryngeal Mask Airway features a cushioned inflatable cuff with reinforced wire coil to provide a direct, sealed route to the lungs. The Laryngeal Mask Airway is ideal for use by anesthesiologists administering general anesthesia during surgery and other medical procedures. The LMA also provides a simpler alternative to endotracheal intubation in emergency situations, allowing for relatively easy insertion and blind intubation, even by novice operators.

- Inflatable cuff fills oral pharynx to help protect airway
- Soft materials and contoured design help reduce throat irritation
- Ideal for use during the administration of anesthesia
- Allows for relatively easy placement and blind intubation
- Translucent tubing and depth marks allow for accurate placement
- Reinforced with wire coil for optimal support and reduced risk of kinking
- Single-use mask is sterile and x-ray detectable

The **Reinforced** masks feature a **wire coil** for maximum strength, reducing risk of kinking



Inflatable Cuff provides direct sealed routed to lungs



<p>1.0 mm Reinforced</p> <p>REORDER NO. 36285</p> <p>NEONATES/INFANTS 5KG</p> <p>PACKAGE QUANTITY 5/Case</p>	<p>1.5 mm Reinforced</p> <p>REORDER NO. 36286</p> <p>INFANTS 5-10KG</p> <p>PACKAGE QUANTITY 5/Case</p>	<p>2.0 mm Reinforced</p> <p>REORDER NO. 36288</p> <p>INFANTS/CHILDREN 10-20KG</p> <p>PACKAGE QUANTITY 5/Case</p>	<p>2.5 mm Reinforced</p> <p>REORDER NO. 36289</p> <p>CHILDREN 20-30KG</p> <p>PACKAGE QUANTITY 5/Case</p>	<p>3.0 mm Reinforced</p> <p>REORDER NO. 36290</p> <p>CHILDREN 30-50KG</p> <p>PACKAGE QUANTITY 5/Case</p>	<p>4.0 mm Reinforced</p> <p>REORDER NO. 36291</p> <p>ADULTS 50-70KG</p> <p>PACKAGE QUANTITY 5/Case</p>	<p>5.0 mm Reinforced</p> <p>REORDER NO. 36292</p> <p>ADULTS 70-100KG</p> <p>PACKAGE QUANTITY 5/Case</p>
---	---	---	---	---	---	--

<p>1.0 mm Non-Reinforced</p> <p>REORDER NO. 36293</p> <p>NEONATES/INFANTS 5KG</p> <p>PACKAGE QUANTITY 10/Case</p>	<p>1.5 mm Non-Reinforced</p> <p>REORDER NO. 36294</p> <p>INFANTS 5-10KG</p> <p>PACKAGE QUANTITY 10/Case</p>	<p>2.0 mm Non-Reinforced</p> <p>REORDER NO. 36295</p> <p>INFANTS/CHILDREN 10-20KG</p> <p>PACKAGE QUANTITY 10/Case</p>	<p>2.5 mm Non-Reinforced</p> <p>REORDER NO. 36296</p> <p>CHILDREN 20-30KG</p> <p>PACKAGE QUANTITY 10/Case</p>	<p>3.0 mm Non-Reinforced</p> <p>REORDER NO. 36297</p> <p>CHILDREN 30-50KG</p> <p>PACKAGE QUANTITY 10/Case</p>	<p>4.0 mm Non-Reinforced</p> <p>REORDER NO. 36298</p> <p>ADULTS 50-70KG</p> <p>PACKAGE QUANTITY 10/Case</p>	<p>5.0 mm Non-Reinforced</p> <p>REORDER NO. 36299</p> <p>ADULTS 70-100KG</p> <p>PACKAGE QUANTITY 10/Case</p>
--	--	--	--	--	--	---



1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2022 Dynarex Corporation

R220622