



# NASOPHARYNGEAL AIRWAYS & KITS

Dynarex Resp-O<sub>2</sub>™ Nasopharyngeal Airways are designed for nasal insertion to help create and maintain an open and secure airway. The Nasopharyngeal Airway features a fixed flange end for safer insertion and internal striations for optimal suctioning efficiency. The Nasopharyngeal Airway Kits include an assortment of airways in multiple sizes for optimal versatility and packets of lubrication jelly for added convenience.

- Designed to be inserted into nasal passage
- Fixed flange design limits the depth of insertion for safe, secure placement
- Constructed from soft materials for enhanced comfort
- Features rounded tip for gentle insertion
- Strong, flexible tubing resists collapse and kinking
- Not made from natural rubber latex



**5.0 mm**

Nasopharyngeal Airway

**REORDER NO. 36275**

20 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**5.5 mm**

Nasopharyngeal Airway

**REORDER NO. 36276**

22 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**6.0 mm**

Nasopharyngeal Airway

**REORDER NO. 36277**

24 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**6.5 mm**

Nasopharyngeal Airway

**REORDER NO. 36278**

26 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**7.0 mm**

Nasopharyngeal Airway

**REORDER NO. 36279**

28 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**7.5 mm**

Nasopharyngeal Airway

**REORDER NO. 36280**

30 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**8.0 mm**

Nasopharyngeal Airway

**REORDER NO. 36281**

32 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**8.5 mm**

Nasopharyngeal Airway

**REORDER NO. 36282**

34 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

**9.0 mm**

Nasopharyngeal Airway

**REORDER NO. 36283**

36 FR  
 STERILE  
 PACKAGE QUANTITY  
 10/Box

Kit of 6

**REORDER NO. 36270**

AIRWAY SIZES INCLUDED (MM)  
 5.5, 6.0, 6.5, 7.0, 7.5, 8.0

6 DYNALUBE PACKETS

PACKAGE QUANTITY  
 1 Kit/Box

Kit of 9

**REORDER NO. 36271**

AIRWAY SIZES INCLUDED (MM)  
 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0

9 DYNALUBE PACKETS

PACKAGE QUANTITY  
 1 Kit/Box



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2022 Dynarex Corporation

R220803