



## SURGICAL SUTURES


Designed to efficiently hold bodily tissues together after surgery or injury, these premium sutures feature high-tensile strength for optimal reliability.

- Premium sutures offer high-tensile strength for optimal reliability
- Designed to provide smooth passage through tissue
- Offers maximum flexibility and consistent knot strength
- Sterile
- 12 / box



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |



NEEDLE KEY

SYM.	TYPE
▽	Precision Reverse Cut 3/8 circle
▼	Reverse Cut 3/8 circle

SURGICAL SUTURES

NYLON SUTURES - SYNTHETIC

18" BLACK MONOFILAMENT | NON-ABSORBABLE

ITEM #	DESC.	CODE	SYM.	SIZE
9100	4-0	PC31	▽	19 mm
9101	6-0	C3	▽	13 mm
9102	5-0	C6	▼	19 mm
9103	4-0	C6	▼	19 mm
9104	3-0	C7	▼	24 mm
9106	5-0	C3	▽	13 mm
9107	4-0	C3	▽	13 mm



POLYPROPYLENE SUTURES - SYNTHETIC

18" BLUE MONOFILAMENT | NON-ABSORBABLE

ITEM #	DESC.	CODE	SYM.	SIZE
9119	5-0	C3	▽	13 mm



POLYGLYCOLIC ACID (PGA) SUTURES

COATED SYNTHETIC

VIOLET | 18" BRAIDED | ABSORBABLE

ITEM #	DESC.	CODE	SYM.	SIZE
9130 *	4-0	C6	▼	19 mm
9131	5-0	C3	▽	13 mm
9132	4-0	C3	▽	13 mm

\* ITEM # 9130 is 30"



SILK SUTURES - NATURAL MATERIAL

BLACK | 18" BRAIDED | NON-ABSORBABLE

ITEM #	DESC.	CODE	SYM.	SIZE
9142	4-0	C6	▼	19 mm
9143	3-0	C6	▼	19 mm



NEEDLE COMPATIBILITY

SYM.	SIZE	DYNAREX	CP MEDICAL	ETHICON	COVIDIEN	LOOK
▽	19 mm	PC31	NPS-2	PS-2	P-12	PC31
▽	13 mm	C3	NP-3	P-3	P-13	C3
▼	19 mm	C6	NFS-2	FS-2	C-13	C6
▼	24 mm	C7	NFS-1	FS-1	C-14	C7
▽	11 mm	C1	NP-1	P-1	P-10	C1

