



WE CARE LIKE FAMILY™

Tillotson True Advantage® Nitrile Exam Gloves



Discover the True Advantage

Manufactured using a patented non-sulfur based curing process which connects the polymer molecules unlike the curing methods of other nitrile gloves. Allows the polymer more freedom of movement, allowing the glove to conform to the hand rather than constrict it.

- The first incredibly soft nitrile glove that improves upon the feel of latex.
- Non-Latex Powder Free Exam Glove with chlorinated lining.
- Meets NFPA standards and is chemotherapy approved.



Reorder No.	Description	Quantity
6501	Nitrile Exam Glove - XS	10/100/Cs
6502	Nitrile Exam Glove - S	10/100/Cs
6503	Nitrile Exam Glove - M	10/100/Cs
6504	Nitrile Exam Glove - L	10/100/Cs
6505	Nitrile Exam Glove - XL	10/100/Cs

True Advantage glove – Quick Facts			
Color	Thickness	Textured	Tensile Strength
	5 mil	Yes	14 MPa, Min.
Cuff	Meets Standards		Market Position
Beaded	ASTM D6319, ASTM D5151 ASTM D6978-05		KC Premium

Meets NFPA standards Thickness may vary



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com