

Medical Imaging

Assistance

- Electroconductive gel media used to facilitate the ultrasonic pulsed echo imaging system
- Acoustically correct for a broad range of frequencies
- Aqueous, will not stain clothing or damage transducers
- Can be used with any equipment that requires an aqueous, conductive gel (for external use only)
- Fragrance free
- Bacteriostatic and non-irritating
- Non-greasy
- Easy refill
- Not made with natural rubber latex

Reorder No.	Description	Quantity
1241	Ultrasound Gel - 0.25 liter (8.5 fl. oz.) Blue	12/Cs
1243	Ultrasound Gel - 1.3 gal (5 liters) Blue	4/Cs
1245	Ultrasound Gel - 0.25 liter (8.5 fl. oz.) Clear	12/Cs
1247	Ultrasound Gel - 1.3 gal (5 liters) Clear	4/Cs

See our website for our DynaLube (#1250 - 1254) & Stethoscopes (#7110 - 7111)















dynarex[®] is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com