



The Monitoring Station provides real-time patient data monitoring for nursing requirements. Functions like nursing administration and alarm setting and remote NIBP measurement help to simplify tedious work for caregivers. All patients' data are well-organized and available to be transferred to other Monitoring Station, and can be checked on Viewer Station, Mobile Station and Web Observer.



The Viewer Station can be located anywhere as needed, such as doctor office, duty office or even corridor. Even if caregivers are out of the department, Mobile Station can provide an uninterrupted monitoring solution. Without installing any software, Web Observer enables caregivers to focus on the patient's condition on the browser.



With the cross server connectivity, MFM-CMS central monitoring system allows supervisors observe the alarm information from different sectors or departments simultaneously, ensuring the safety of all patients at one glance.



MFM-CMS Central Monitoring System

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



A world of potential

Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com
U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

Cross departments, Cross devices

MFM-CMS Distributed Monitoring System



ENG-PM-MFM-CMS III
V1.0-20230329



A world of potential

Bring the data to your sight

The new generation MFM-CMS brings you a comprehensive bi-directional central monitoring solution for all EDAN monitors. It offers seamless monitoring management from transportation to hospitalization, from mobile patients to bedside cares, from intensive units to operation theatres, from in-hospital to out-of-hospital, bringing easy access to well-organized monitoring information.



Max. 4 screens




Max.128 patients/sector



Data review



Remote control



Storage capacity for each patient:

720 h	720 h	720 h	720 h	720 h
NIBP Review	Alarm Event	ARR Review	Trend	Full-disclosure Waveform

Patient-oriented medical record management

New generation MFM-CMS can integrate the 4 product solutions of Edan and simultaneously view the patient physiological data provided by ventilators, anesthesia machines and infusion pumps of mainstream manufacturers. Patient's reports, including other treatment tools, can all be checked on MFM-CMS. By supporting POC Test reports, Holter analysis reports and ambulatory BP reports, doctors can have an overview of patients' conditions in specific departments, especially the cardiology department, ICU and OR.



Edan:

Blood gas and chemistry analyzer: i15




Magnetic immunoassay system: m16



Other Products:

Infusion & Syringe pump: Medcaptain HP30, HP60



Nursing system for spot check

For both long-time monitoring and spot check application, MFM-CMS can provide caregivers with a practical solution, including Patient List, Nursing Record Sheet, TEMP Record Sheet and Work Shift Report, which relieves nurses from heavy paperwork. With the help of MFM-CMS, users will need no more extra nursing administration system to adapt to the spot check scenario. Easy operation and powerful functions make nursing workflow closed-loop.

	Time	Vital S	Intake	Drug	Basic I	Canal	Perfric	Psychc	Pain	Skin	Body I	Spec	Work	Signa	Ope
<input type="checkbox"/>	7-19-2022 10:18:4	Normal	Normal	/	Close A	Normal	Normal	Normal	Normal	/	Lie Flat	/	Jane	Mary	
<input type="checkbox"/>	7-19-2022 09:23:0	Normal	Normal	/	Normal	Normal	Normal	Close A	Normal	/	Prone F	/	Mary	Emily	


Work shift report:

Release nurses from heavy paperwork and make all work shift reports in good order.

	Ope	Time	TEMP(°C)	NIBP DIA(mmHg)	NIBP SYS(mmHg)	NIBP MAP(mmHg)	PR(bpm)	HR(bpm)
<input type="checkbox"/>		7-19-2022 10:36:03	36.0	77	120	80	60	
<input type="checkbox"/>		7-19-2022 10:36:52	36.0	71	119	85	60	
<input type="checkbox"/>		7-19-2022 10:38:04	36.0	76	120	88	60	

Nursing record sheet:

All nursing records for the same patient can be displayed and printed out in a tabular format for caregivers to observe changes in the patient's condition.



TEMP record sheet:

Carefully designed icons and clear table draw the intuitive display of patient's TEMP data changes.

