

INTRODUCING THE
HSS SERIES SURGICAL STOOL



HSS500-S SADDLE SEAT STYLE



HSS500-WD WATER DROP SEAT STYLE



HSS500-W WEDGE SEAT STYLE

BUILT FOR PERFORMANCE
AND COMFORT

HAUSTED SURGICAL STOOL

IDEAL FOR EYE, ENT AND PLASTIC PROCEDURES



360° Rotating
Urethane Arm Rest
for Comfortable
Arm Positioning

Adjustable Contoured Backrest

**300
POUNDS**
300 lb (136 kg)
Weight Capacity
EVENLY DISTRIBUTED

Dual-density Foam Seat
with **DURA-GEL**[®] Insert
for Enhanced Comfort



180° Arm Rest Support for more
Stable Arm Positioning


Adjustable Height Range
21 - 28 inches (53 - 71 cm)

OPTIONAL ACCESSORY:
Riser Kit Increases Height Range
an additional 2 inches (5 cm)
Item No: H081025

ABS Thermoformed
Hood for Easy Cleaning



Complete the HSS Series Surgical Stool
Experience with the new Hausted
ESC² Ophthalmology Platform
Item No: ESC2EYEST



FOOT-CONTROLLED HYDRAULICS
with Three-Spindle Base and
Four 3 inch (8 cm) Dual Tente[®]
Casters for Greater Stability



Front
Locking
System

ENGINEERED TO PROVIDE
ENHANCED
SURGEON COMFORT



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.



EC REP

CE Partner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.eu

**graham
field**