



A comprehensive platform scale that supports up to 1000 pounds offers “live” handrails and a low profile platform to provide a stable and comfortable weighing experience for patients of any size. Ideal for any type of facility including primary care, rehab or hospital.

Advantages

- ✓ **EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **“Live” Handrails**
Advanced technology allows the patient to hold onto handrails for stability while scale accurately calculates weight
- ✓ **Low Profile Platform**
Low platform height requires only minimal step up, offering patient a comfortable and stable weighing experience
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient’s BMI
- ✓ **180° Swivel Display**
The display head has the ability to pivot for the convenience of the caregiver, allowing line-of-sight viewing and information privacy
- ✓ **On-Screen Help Menu**
Guides users through operating functions and provides easy access to multiple settings
- ✓ **Interactive Display**
High-contrast color screen provides a user-friendly interface
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG LockOut, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost
- ✓ **Everlock® (patent pending)**
Option that allows users to permanently lock the scale in either LB or KG to improve safety



PROPLUSROD:
Optional Height Rod



180° Swivel Display



Specifications

- Capacity: 1000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
22" x 15 3/4" x 2 1/2" / 559 mm x 400 mm x 64 mm
- Product Footprint (w x d x h):
31 1/4" x 27 1/4" x 49 3/4" / 794 mm x 692 mm x 1264 mm
- Product Weight: 52 lb / 23 kg
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Wheels: 2
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Pelstar wireless technology
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Everlock®
- 2 year limited warranty

Model	Description
1100KL	Digital platform scale, pounds and kilograms
1100KL-BT	1100KL with built-in Pelstar wireless technology
1100KG	Digital platform scale, kilograms only
1100KG-BT	1100KG with built-in Pelstar wireless technology
1100KLHR	1100KL with optional height rod
1100KGHR	1100KG with optional height rod
SS-1100KL	ScaleSurance 2 year extended warranty (1100 models)*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PROPLUSROD	Optional height rod
PELSTARONE-P	Wireless thermal printer