



This classic “Bigfoot” scale has a deep platform with non-slip surface and an 8” easy-to-read raised dial—good for patients with vision challenges or those who prefer more foot room.

### Advantages

- ✓ **400 lb / 180 kg Capacity Floor scale**  
The high capacity of this floor scale provides more versatile use
- ✓ **8" Raised Dial Display**  
Large, easy-to-read dial display is raised above the scale base to allow for more footroom on the platform
- ✓ **Weighs in LB/KG**  
Dial display shows weight in both pounds and kilograms providing versatility



**Specifications**

- Capacity: 400 lb /180 kg
- Graduation: 1 lb / 1 kg
- 8" dial (raised 4")
- Platform Size (w x d x h):  
11" x 12 ½" x 3" / 279 mm x 318 mm x 76 mm
- Product Footprint (w x d x h):  
11 ¾" x 18 ¼" x 7 ¼" / 298 mm x 464 mm x 184 mm
- Product Weight: 9 lb / 4 kg
- 1 year limited warranty

Model	Description
160KL	Mechanical floor scale - kilograms & pounds (2 pack)
160LB	Mechanical floor scale - pounds only (2 pack)
160KG	Mechanical floor scale - kilograms only (2 pack)
160KLS	Mechanical floor scale - kilograms & pounds (single)
160LBS	Mechanical floor scale - pounds only (single)
160KGS	Mechanical floor scale - kilograms only (single)