




2210KL



2210KL4

The 2210KL neonatal scale was designed for use in the NICU, nursery, delivery room and pediatric practices where high degree of accuracy is critical to ensure that babies are gaining weight (2210KL4 comes with a four sided tray). A convenient cart is available to easily move the scale from location to location.

Advantages

- ✓ **Motion Detection**
Motion sensing technology compensates for baby's movement and high-resolution readings provide more accurate weight measurements, ideal for the neonatal environment
- ✓ **Value-Added Features**
Operating functions (LB / KG / GM Lockout, LB / GM Conversion, LB / KG Conversion, Zero, Tare, Reweigh / Recall, Auto Zero, Auto Off) included at no additional cost
- ✓ **Convenient Features**
Plastic tray is easy-to-clean and heavy-duty built-in handles allow for portability
- ✓ **Built-In Measuring Tape**
Easy-to-read measuring tape allows user to measure both weight and height simultaneously



2210KL / 2210KL4 Specifications



- Capacity: 45 lb / 20 kg / 20,000 gm
- Resolution: 1 g (to 10,000 g); 5 g (10,001 to 20,000 g); 0.1 oz (0 – 45 lb)
- Display Size (w x h): 5 ¼" x 1 ⅛" / 133 mm x 35 mm, LCD display
- Functions: LB / KG / GM Lockout, LB / GM Conversion, LB / KG Conversion, Tare, Zero, Reweigh / Recall, Auto Zero, Auto Off, Motion Detection
- Removable curved clear acrylic tray
- Tray Size (2210KL) (w x d x h): 24" x 15 ⅜" x 3" / 610 mm x 391 mm x 76 mm
- Tray Size (2210KL4) (w x d x h): 27" x 15 ⅜" x 3" / 686 mm x 391 mm x 76 mm
- Product Weight (2210KL): 24 lb / 10.8 kg
- Product Weight (2210KL4): 25 lb / 11.3 kg
- Measuring Tape: 0" - 23" / 0 cm - 60 cm
- Heavy duty handles for ease of moving
- 6 C batteries included
- 2 year limited warranty - eligible for ScaleSurance extended warranty

2210CART Specifications

- Capacity: 200 lb / 90 kg
- Top shelf has handle for easy grip. Non-conductive wheels (2 wheels with locks)
- Durable, easy-to-clean storage shelf
- Cart Dimensions: 26 ½" (w) x 18" (d) x 34 ⅛" (h) / 673 mm (w) x 457 mm (d) x 867 mm (h)
- Lockable Cabinet: 23" (w) x 16 ¾" (d) x 17" (h) / 584 mm (w) x 427 mm (d) x 432 mm (h)
- Use with neonatal / pediatric scales: 2210KL, 2210KL4, 553KL, 549KL
- 1 Year Limited Warranty



Model	Description
2210KL	High resolution digital neonatal / pediatric tray scale
2210KL4	High resolution digital neonatal / pediatric 4 sided tray scale
2210KLCT	2210KL with rolling cart included
2210KL4CT	2210KL4 with rolling cart included
SS-2210KL	ScaleSurance 2 year extended warranty for 2210KL/2210KL4*
2210CART	Rolling cart for all models above
BLUETOOTHONE-P	Wireless thermal printer