



The deep infant seat on the 524KL reduces the movement of young children to speed up the weighing process. Weight is captured easily while the seat supports the child's head and body.

Advantages

- EMR Connectivity**
 Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- 1" / 25 mm LCD Display**
 The large display shows the user an easy-to-read weight measurement
- Motion Detection**
 Scale has the ability to weigh accurately despite movement of the subject on the scale
- Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- Convenient Weighing Seat**
 Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement



Specifications

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb 0.2 oz / 20 - 50 lb 0.5 oz 0 - 9 kg / 5g;
9 - 23 kg / 10 g
- Connectivity via USB
- 1" / 25 mm LCD display
- Functions: LB/KG Conversion, LB/KG Lock Out, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- Seat Size (w x d x h):
24 1/4" x 14 1/4" x 18 1/2"
616 mm x 362 mm x 470 mm
- Product Footprint (w x d x h):
24 1/2" x 21" x 23 1/8"
622 mm x 533 mm x 587 mm
- Product Weight: 21 lb / 10 kg
- 6 AA batteries included, Optional 120V adapter (order # ADPT31)
- 2 year limited warranty

Model	Description
524KL	Digital pediatric seat scale, pounds & kilograms
SS-524KL	ScaleSurance 2 year extended warranty
ADPT31	Optional 120V adapter
BLUETOOTHONE-P	Wireless thermal printer