



This scale is ideal for infants and with its digital height rod practitioners can measure both weight and height simultaneously.

Advantages

- ✓ **EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Motion Detection**
Scale has the ability to weigh accurately despite movement of the subject on the scale
- ✓ **Scale Value-Added Features**
Operating functions (LB/KG Conversion, KG Lock Out, Zero Out/Tare, Hold/Release, Auto Off) included at no additional cost
- ✓ **Digital Height Rod Included**
Easy to read height rod allows user to measure both weight and height simultaneously



Scale Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h):
24 7/8" x 13" x 1 3/8"
632 mm x 330 mm x 35 mm
- Display: 1" / 25 mm LCD
- Power Source: 6 AA batteries included,
optional power adapter
(order # ADPT30)
- Connectivity: USB, optional Pelstar
wireless technology
- Functions: Zero, Tare, Hold / Release,
Auto Off
- 2 year limited warranty

Height Rod Specifications

- Measurement Range:
13 3/4" – 31 1/2" / 35 -80 cm
- Resolution: 1/16" / 0.1 cm
- Power Source: 2 AAA batteries (included)
- Display: LCD, 5/8" / 16 mm digit height
- Digital display provides height in
centimeters, inches, & feet
- Functions: Auto Off, Hold
- 2 year limited warranty

| Model | Description |
|--------------|--|
| 553KL-EHR | Digital pediatric scale (pounds & kilograms) with digital height rod |
| 553KG-EHR | Digital pediatric scale (kilograms only) with digital height rod |
| 2210CART | Rolling cart for 553KL-EHR & 553KG-EHR |
| ADPT30 | Optional 100-240V adapter for scale |
| PLUGSET30 | Plug adapters for use outside North America for scale |
| PELSTARONE-P | Wireless thermal printer for scale |