



The 553KG is perfect for newborns, infants and toddlers. The scale operates on batteries or AC power so you can use it anywhere. Practices benefit from one-button operation that makes it easy to use while ensuring that the baby remains safely on the scale.

### Advantages

- EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Value-Added Features**  
Operating functions (Zero, Tare, Hold/Release, Auto Off) included at no additional cost
- Motion Detection**  
Scale has the ability to weigh accurately despite movement of the subject on the scale
- Measuring Tape Included**  
Easy-to-read measuring tape allows user to measure both weight and height simultaneously



Carrying Case  
Item # 553CASE



**Specifications**

- Capacity: 20 kg
- Resolution: 10 g
- Tray Size (w x d x h):  
24 7/8" x 13" x 1 3/8"  
632 mm x 330 mm x 35 mm
- Product Footprint (w x d x h):  
24 7/8" x 16 1/2" x 3 7/8"  
632 mm x 419 mm x 98 mm
- Product Weight: 10 lb / 5 kg
- Measuring Tape: 0" - 24" / 0 cm - 58 cm
- Display: 1" / 25 mm LCD
- Power Source: 6 AA batteries included, optional power adapter (order # ADPT30)
- Connectivity: USB, optional Pelstar wireless technology
- Carrying Case: Optional carrying case (order # 553CASE)
- Functions: Zero, Tare, Hold / Release, Auto Off
- 2 Year Limited Warranty

Model	Description
553KG	Digital pediatric tray scale, kilograms only
553KG-HR	553KG with attachable mechanical height rod
553KG-EHR	553KG with attachable digital height rod
SS-553KL	ScaleSurance 2 year extended warranty for 553KG
553HR	Attachable mechanical height rod for 553KG
553EHR	Attachable digital height rod for 553KG
2210CART	Rolling cart for 553KG
553CASE	Carrying case for 553KG
ADPT30	Optional 100-240V adapter
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer