



Advantages

- ✓ **400 lb/180 kg Capacity Floor Scale**
The high capacity of this floor scale provides more versatile use
- ✓ **1 ½" LCD Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **5-Second Weight Hold**
Weight measurement will stay displayed for 5 seconds after stepping off scale allowing the user time to stabilize themselves or the patient

Health o meter®

P r o f e s s i o n a l



Carrying Case for 800KL
Item# 64771



Specifications

- Capacity: 400 lb / 180 kg
- Resolution: 0.1 lb / 0.05 kg
- 1 ½" LCD display
- Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off
- Platform Size (w x d x h):
13 ¾" x 10 ¾" x 1 ½"
349 mm x 273 mm x 38 mm
- Product Footprint (w x d x h):
14 ¼" x 11 ⅞" x 1 ½"
362 mm x 302 mm x 38 mm
- (2) 3-Volt Lithium battery included
- Optional Carrying Case Available (order # 64771)
- 1 Year Limited Warranty