

Health o meter®

P r o f e s s i o n a l



Advantages

- 440 lb / 200 kg Capacity Floor Scale**
The high capacity of this floor scale provides more versatile use
- 1-1/4" LCD Display**
The large display shows the user an easy-to-read weight measurement
- EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience

Health o meter[®]
P r o f e s s i o n a l



844KL Display



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
844KL	2	15.125" x 14.75" x 6" / 0.77	8.5 lb	10892076002401



Specifications

- Capacity: 440 lb / 200 kg / 31 st
- Resolution: 0.1 lb / 0.05 kg / 1/8 lb
- Platform Size: 12 5/8" (w) x 12 5/8" (d) x 1 7/8" (h) / 321 mm (w) x 321 mm (d) x 48 mm (h)
- Display: 1 1/4" LCD
- Power Source: (2) 3-Volt Lithium batteries included
- Connectivity: USB
- Optional Carrying Case item # 64771
- Functions: LB / KG / ST Switch, Auto Zero,
- 5 Second Weight Hold, Auto Off
- 1 Year Limited Warranty



Carrying Case for 844KL
Item# 64771

Health o meter[®]
P r o f e s s i o n a l
www.homscales.com

Health o meter[®] Professional Scales
9500 West 55th Street USA
McCook, IL 60525-7110
Phone: 1-800-815-6615
Fax: 1-800-654-7330