



High Capacity:
894KLT: 400lb / 180kg
895KLT: 550lb / 250kg

Advantages

- ✓ **Bilingual Audible Weight Measurement**
For the visually impaired, this “talking” scale announces the weight measurement in English or Spanish
- ✓ **20-Second Weight Hold**
Weight measurement will stay displayed for 20 seconds after stepping off scale allowing the user time to stabilize themselves or the patient
- ✓ **1 ½" LCD Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **No Carpet Feet Needed**
Outperforms all leading brands during use on carpet.

Health o meter®

P r o f e s s i o n a l



Optional
Carrying Case
Item# 64771



Specifications

- 894KLT Capacity: 400 lb / 180 kg
- 895KLT Capacity: 550 lb / 250 kg
- Resolution: 0.2 lb / 0.1 kg
- 1 ½" LCD display
- Functions: LB / KG Switch, Auto Zero, Auto Off, 20 Second Weight Hold, Talking Scale, English / Spanish Button, Voice Disable Option
- Platform Size (w x d x h):
13 ¾" x 10 ¾" x 1 ½"
349 mm x 273 mm x 38 mm
- Product Footprint (w x d x h):
14 ¼" x 11 ⅞" x 1 ½"
362 mm x 302 mm x 38 mm
- Power Source: 2 AAA batteries included
- Optional carrying case available (order # 64771)
- 1 year limited warranty

Health o meter®
Professional
www.homscales.com