



InteliSwab™ COVID-19 Rapid Test makes testing easy.

A remarkably simple COVID-19 rapid antigen test

- Self-collected, visually-read lateral flow test
- Both symptomatic and asymptomatic¹ use
- Intuitive workflow: just “swab, swirl, and see”
- Gentle lower nasal swab
- < 1 minute of “hands-on” time, 30 min time to result

3 InteliSwab™ configurations:

- Over-the-counter Home Test (2-pack serial test)
- Prescription Home Test (1-pack)
- Point-of-care bulk packs (healthcare professional use)

Ages 15+

NOW FDA AUTHORIZED

¹Asymptomatic use for InteliSwab™ COVID-19 Rapid Test and Rapid Test Pro only



Dramatically simplify COVID-19 testing with IntelliSwab™

Simple and easy-to-use

- Just 3 simple steps: swab, swirl and see your result
- Intuitive workflow has no difficult or confusing steps
 - ✓ No dropper and precise number of drops required
 - ✓ No complex manipulation of swab
 - ✓ No alignment of test strip onto IFU for interpretation
 - ✓ No risk of damaging unprotected assay strip with accidental touching
- No batteries, instruments or smartphones required to get a result
- Clinically proven easy-to-use

Quick and efficient

- Workflow requires < 1 minute of active, hands-on time
- IntelliSwab™ provides test results in hand in just 30 minutes
- No lab send out needed


Practical

- 9-month shelf life on all tests
- No need for ancillary supplies such as nasal swabs or readers

Accurate

- Overall test performance: Positive Percent Agreement 84%, Negative Percent Agreement 98%
- Clinical performance proven with untrained and unproctored lay-users

Testing will be critical throughout all stages of epidemic transition



Configuration	IntelliSwab™ COVID-19 Rapid Test	IntelliSwab™ COVID-19 Rapid Test Rx	IntelliSwab™ COVID-19 Rapid Test Pro	
	Over-the-Counter	Prescription Home Test	Point-of-Care	
Catalog Number	1001-0616	1001-0620	1001-0614	1001-0615
Test devices	24 x 2-test kits	24 x 1-test kits	25 tests	100 tests