

GUARDIANPLUS

Sedation Monitor

The Guardian Plus monitor has an exceptional user-interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Experience the immediate and measurable benefits of the smartphone-like user-interface that delivers 90% of all common monitoring functions with just two steps.



7 LBS

Weight

4 HR

Battery
Life

8

Waveforms

10 IN

Touch Screen
Display

Features

- 10" touchscreen with gesture control
- Customizable alarms and thresholds
- Adult, pediatric, and neonate presets
- 18 configurable short cuts and quickkeys
- Integrated Co2 option

Standard Parameters:

ECG (3-5 Lead), HR, RESP, Masimo SET SpO2, NIBP, PR, TEMP

Optional Parameters:

Nellcor OxiMax SPO2, CO2

Clinician-focused Technology for Enhanced Patient Care

- Parameter list tile with automatic highlight of abnormal readings
- One touch to drag icon tab right and instantly view minitrend screen or set baseline vitals at the start of a patient case
- Remote view and alarm watch bed-to-bed allows up to 12 patients to be viewed
- Serial DIAP and HL7 communication protocols provide seamless transfer of data to EMR



Primary Care



ASC





Easy Transport



Pole Mount



Bedside Mount

Ergonomic Design

Delivers Exceptional Performance and Reliability:

- Manufactured with industry-leading plastics that stand up to the most stringent medical grade cleaners
- IP44 rated, resistant to damage often caused by continuous cleaning and usage
- The concealed grab-and-go handle eliminates hard to clean crevices and eases clinical workflow while on the move
- Flexible mounting for bedside and transport provides options to fit varying care environments

Unmatched Visual Experience

- Smartphone-like multi-gesture, capacitive touchscreen display
- HD display with wide angle view $\geq 170^\circ$
- Flat user interface design, common across platform
- Reach 90% of common monitor functions in 2 steps

Accessories

Accessories include one battery, one adult cuff, 3.5 m NIBP hose, one line cord, one roll of paper, ECG N/T cable, 3-lead adult/pediatric snap lead wire set, Masimo or Nellcor SpO2 cable and sensor, temperature probe adapter cable, quick reference guide, and one operator's manual. **All other accessories must be ordered separately.**

BATTERY & CABLES

- Lithium Ion Battery — **115-065140-00**
- Power Cord, 110 Volt — **0012-25-0001**
- Serial Port Cable — **009-009117-00**

ECG ACCESSORIES

- ECG Trunk Cable, 10' (3.1 m), Defibrillation-proof, for N/T — **009-005266-00**
- ECG Trunk Cable, 20' (6.2 m), Defibrillation-proof, for N/T — **009-005267-00**
- ECG Trunk Cable, 10' (3.1 m), ESU-Proof, for N/T — **009-005268-00**
- ECG Trunk Cable, 20' (6.2 m), ESU-Proof, for N/T — **009-005269-00**
- 3-Lead, N/T, AHA, pinch/clip, 24" — **009-004771-00**
- 3-Lead, N/T, AHA, pinch/clip, 36" — **009-004772-00**
- 3-Lead, N/T, AHA, snap, 24" — **009-004765-00**
- 3-Lead, N/T, AHA, snap, 36" — **009-004766-00**
- 5-Lead, N/T, AHA, pinch/clip, 24" — **009-004786-00**
- 5-Lead, N/T, AHA, pinch/clip, 36" — **009-004787-00**
- 5-Lead, N/T, AHA, snap, 24" — **009-004782-00**
- 5-Lead, N/T, AHA, snap, 36" — **009-004783-00**

PAPER

- Paper, Thermal Strip Local Printer (12/rolls) — **0683-00-0505-02**

SPO2 ACCESSORIES - MASIMO RD SET

Reusable

- SpO2 sensor, RD SET DCI, Adult (> 30 kg) — **040-003376-00**
- SpO2 sensor, RD SET DCI, Pediatric (10-50 kg) — **040-003377-00**
- SpO2 sensor, RD Set TC-I, tip-clip ear sensor, 3ft, Adult (> 30 kg) — **040-003380-00**

Disposable

- SpO2 sensor, RD SET Adt, adhesive, Adult (> 30 kg, pkg/20) — **040-003382-00**
- SpO2 sensor, RD SET PDT, adhesive, Pediatric (10 to 50 kg, pkg/20) — **040-003383-00**
- SpO2 sensor, RD Set Inf, adhesive, Infant (3 to 20 kg, pkg/20) — **040-003384-00**
- SpO2 sensor, RD Set Neo, adhesive, Neonatal /Adult (<3kg or >40kg, pkg/20) — **040-003385-00**
- SpO2 sensor, RD Set NeoPt, adhesive, Neonatal (<1 kg pkg/20) — **040-003386-00**
- SpO2 sensor, RD Set NeoPt-500, non-adhesive, (<1 kg pkg/20) — **040-003387-00**

Cables

- Masimo RD SET cable, 8-pin — **040-003310-00**
- Masimo RD to LNCS adapter cable — **040-003381-00**
- Masimo LNCS to RD adapter — **040-003426-00**

SPO2 ACCESSORIES, NELLCOR

- Nellcor SpO2 extension cable, 8pin — **0010-20-42712**

Accessories (Continued)

NIBP ACCESSORIES

- Nibp Tubing, Adu/Ped/Neo, with Connectors (3m) — **6200-30-09688**
- Child, reusable cuff, 10-19cm (limb) — **115-027713-00**
- Small Adult, reusable cuff, 18-26cm (limb) — **115-027714-00**
- Adult, reusable cuff, 25-35cm (limb) — **115-027715-00**
- Adult, long reusable cuff, 25-35cm (limb) — **115-027718-00**
- Large Adult, reusable cuff, 33-47cm (limb) — **115-027716-00**
- Large Adult, long reusable cuff, 33-47cm (limb) — **115-027719-00**
- Adult Thigh, reusable cuff, 46-66cm (limb) — **115-027717-00**
- Reusable cuff starter kit — **115-031807-00**

MOUNTING

- Wall Mount with Quick Release - M Series — **045-003427-00**
- Rolling Stand with Quick Release Mount, Cable Management, One Storage Basket — **045-002960-00**
- Storage Basket, Additional (2nd) Without Handle — **045-002936-00**
- Cleaning/disinfecting Wipes Mount — **045-002935-00**
- Bedrail Hook — **8000-30-90169**
- Bedrail Clamp — **8000-30-90170**
- Cable Hooks / Attach to M-series Or Vhm Arm — **042-007382-00**
- Quick Release Mount for Arm Mount and Rolling Stand — **045-003424-00**

CO2 ACCESSORIES

- Dryline II Water Trap, Adult/Pediatric (10/box) — **115-058733-00**
- Dryline II Water Trap, Neonatal (10/box) — **115-058734-00**
- Dryline Airway Adapter, Adult/Pediatric, Straight (10/box) — **115-043020-00**