Traction / Decompression Systems and Therapeutic Tables

Neural-Flex® Decompression System—ME4100



Features

Lateral motion decompression drive mount

Powered front and rear sections smoothly transition the patient to a centralized position

Adjustable tension spring

Bolster-less leg elevation

Powered flexion and extension to 45 degrees

Comfortable headrest with large breathe-through opening

Full body vibration

Retractable arm rests allow more

surface for larger patients

Upper section pivots at L-5

Integrated powered caudal surfaces and elevation control 18 in to 39 in (46 cm to 99 cm)

. . . .

Locking casters



Table Specifications

 $110 \ V{\sim}\ 50/60 \ Hz \ or \ 220 \ V{\sim}\ 50/60 \ Hz$ Dimensions: 98 in (L) x 28 in (W) x 43 in (H) 250 cm (L) x 71 cm (W) x 109 cm (H)

Weight: 251 lbs. (114 kg)

Warranty: 2 years, no warranty on fabric

tears

Included Components:

Neural-Flex table

MTD 4000 programmable traction head

unit

Quick-Ease™ Lumbar belt System

Max Relax Face Cushion

Clinical Cervical Device with RB Pillow

In-Service Videos and Online Training

Colors





Traction / Decompression

MTD 4000, Traction Head Unit—ME4000



Features

Maximum Output Force: 90 Kg / 198 lbs.

Easy to use

Active displays show all treatment parameters

and progress.

Multiple sensors and safety controls

Adjustable Hold/Rest times

Continuous and intermittent traction with

multiple speed selection

Mounts to almost any table platform

Memory to store up to 8 user-defined protocols

Specifications

AC 110 V~120 V, 50/60 Hz or

AC 220 $V\sim$ 240 V, 50/60 Hz

L 14.2 in \times W 12.2 in \times T 9 in

(L 36 cm \times W 31 cm \times T 22.9 cm)

32 lbs. (14.5 kg)

Warranty: 2 years

Mettler Traction System—ME4200



Table Specifications

100~240 VAC, 50/60 Hz

Height: 18.5 in to 35.8 in

(47 cm to 91 cm)

Weight: 227.7 lbs. (103.5 kg)

Warranty: 1 year on Mechanical Parts

2 years on Electronics

3 years on Frame

Features

1 Motorized Actuator for Hi/Lo Adj.

Equipped with 1 Motorized Actuator for Center Arching

3 Sides Elevating Bars

Twin "Oversized" Pillars for Added Stability Individual "Leveling Foot" on All 4 Corners

Retractable Casters at Both Ends

Twin Gas Struts for Solid Head Section

Support

Included Components:


Mettler Traction Table


MTD 4000 programmable traction head

unit

Quick-Ease™ Lumbar belt System

2 x Bolsters







Specialty Therapeutic Tables

7 Section Therapeutic Table—ME4700



100~240 VAC, 50/60 Hz Height: 22 in to 40 in (56 cm to 99 cm)

Weight: 214.5 lbs. (97.5 kg)

Warranty: 1year on Mechanical Parts

2 years on Electronics 3 years on Frame


Colors



Features

1 Motorized Actuator for Hi/Lo Adj.
Twin "Oversized" Pillars for Added Stability
4 Sides Elevating Bars
Individual "Leveling Foot" on All 4 Corners
Retractable Casters at Both Ends
Twin Gas Struts for Solid Leg Section Support

Single Gas Strut for Center Arching & Head Section



8 Section Chiropractic Table—ME4800



Specifications

100~240 VAC, 50/60 Hz Height: 21.3 in to 38.6 in

(54 cm to 98 cm)

Weight: 205.7 lbs. (93.5 kg)

Warranty: 1year on Mechanical Parts

2 years on Electronics 3 years on Frame

Colors



Features

- 1 Motorized Actuator for Hi/Lo Adj.
- 1 Motorized Actuator for Hi/Lo Adj. (Head Section)
- 1 Motorized Actuator for Tilting Adj. (Head Section)
- 1 Joystick for Hi/Lo & Tilting Adj.
- 2 Sides Elevating Bars

Adjustable Tension Drop - Cervical, Dorsal (Thoracic), Lumbar, Pelvic


Twin "Oversized" Pillars for Added Stability

2 Individual Leveling Foot at Head Section

2 Hauling Wheels at the Rear Section

Designed to be lifted and moved around by lifting JUST the Arm Rest

2 Fixed Arm Rests for absolute comfort & reliability



Therapeutic Tables

2 Section Therapeutic Table—ME4500



Specifications

100~240 VAC, 50/60 Hz Height: 18.1 in to 35.8 in (46 cm to 91 cm)

Weight: 150.7 lbs. (68.5 kg)

Warranty: 1year on Mechanical Parts

2 years on Electronics 3 years on Frame

Features

Motorized Actuator for Hi/Lo Adj.


Locking casters

4 Sides Elevating Bars

Individual "Leveling Foot" on All 4 Corners

Twin Gas Struts for Solid Head Section

Support



3 Section Therapeutic Table, no drop end—ME4400



Specifications

100~240 VAC, 50/60 Hz Height: 18.1 in to 35.8 in (46 cm to 91 cm)

Weight: 172.7 lbs. (78.5 kg)

Warranty: 1year on Mechanical Parts

2 years on Electronics3 years on Frame

Features

Motorized Actuator for Hi/Lo Adj.


Locking casters

4 Sides Elevating Bars

Individual "Leveling Foot" on All 4 Corners

Twin Gas Struts for Solid Head Section

Support



3 Section Therapeutic Table—ME4600



Specifications

100~240 VAC, 50/60 Hz Height: 18.1 in to 35.8 in

(46 cm to 91 cm)

Weight: 159.5 lbs. (72.5 kg)

Warranty: 1year on Mechanical Parts

2 years on Electronics

3 years on Frame

Features

Motorized Actuator for Hi/Lo Adj.


Locking casters

Hand control/Foot control 2 Buttons

Individual "Leveling Foot" on All 4 Corners

Twin Gas Struts for Solid Head Section

Support





1333 South Claudina Street • Anaheim, CA 92805 • USA 1.800.854.9305 • 1.714.533.2221 • Fax: 1.714.635.7539 www.mettlerelectronics.com • mail@mettlerelectronics.com