



# CFPM200 & CFPM201 Imaging Tables

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.



## **Motions**

- Motion 1 Height Range 25" 42" (CFPM201)
- Motion 1 Height Range 25" 43" (CFPM200)
- Motion 2 Lateral Tilt ±15° (CFPM200)
- Motion 2 Trendelenburg Tilt +15/-11°
   (CFPM201)

## **Features**

- Choice of integrated headrest or rectangular carbon fiber tops with 56" of metal free imaging
- 2" Radiolucent Comfort Foam® table top pad
- Warranty is 4 years parts and 3 years labor

# **SPECIFICATIONS**

- Large Adjustable Height Range 25"-43" (CFPM200), 25"-42" (CFPM201)
- 56" Metal Free Imaging Space
- Cost Effective 2 Movement Table Models
- 500 lb Patient Weight Capacity
- Warranty: 4 years parts, 3 years labor

CFPM200 #84522 Integrated Rectangular

#84526 Integrated CFPM201 #84528 Rectangular

CFF1V12UU #84524 F	Rectangular CIFIVIZOI #84528 Rectangular	
INCLUDED FEATURES		
Carbon Fiber Top 56" Metal free imaging space	Integrated Headrest 22" (56cm) wide x 84" (214cm) long Rectangular 24" (61cm) wide x 80" (203cm) long	
CFPM200 MOVEMENTS		
Motion 1	Electric Height Range 25" - 43" (64cm-109cm) (without pad)	
Motion 2	Electric Lateral tilt ±15°	
CFPM201 MOVEMENTS		
Motion 1	Electric Height Range 25" - 42" (64-107cm) (without pad)	
Motion 2	Electric Trendelenburg tilt +15°/-11°	
Controls	All functions hand control	
Weight Capacity	500 lb (227kg)	
Product Weight	475 lb (216kg)	
Included Features	2" (5cm) Radiolucent Comfort Foam™ table top pad, 4" (10cm) locking casters, Crescent face rest pad (included with integrated headrest top), patient restraint straps (2)	
Warranty	4 years parts, 3 years labor	
	4 years parts, 5 years labor	
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate	
,	FDA listed, CE marked, cETLus listed, CB Certificate	
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate	
Safety Listings  OPTIONS must be ordered to	FDA listed, CE marked, cETLus listed, CB Certificate  with table.	
Safety Listings  OPTIONS must be ordered a 5"(13cm) casters	FDA listed, CE marked, cETLus listed, CB Certificate  vith table.  Recommended for mobile service providers	

ACCESSORIES order any time.		
Arm Hammock	Use with Integrated Headrest Top	48751
Carbon Fiber Arm Board (6" wide)	180° Articulation	83575
Fluoro Extender	Polycarbonate	33975
Interventional Supine Pillow System	Contoured pad and support wedge pad	81821
Spine Positioning System II	8 piece system	81819
T-rail Clamp	6" (15cm) long rail	63336



CFPM200 Rectangular Top Table

# **ACCESSORIES**

#### SPINE POSITIONING SYSTEM II

Radiolucency & adjustability for ideal cervical, thoracic & lumbar visualization. Minimize unwanted movement & reduce procedure time.



#### ARM HAMMOCK

This unique product will comfortably & easily support your patient's arms when prone. Easily attached to the Integrated Head Rest table top.



#### SINGLE CARBON FIBER ARM BOARD

Fully radiolucent with a 180° range of motion, the arm board requires no mounting or attachment hardware and simply slides under the table top pad. It can be used on any fluoroscopy table & will support arms and feet in a variety of positioning configurations.



#### FLUORO EXTENDER

Made of radiolucent, polycarbonate material, it's organic shape is conducive to a variety of table placement applications.

30" long x 23.5" wide x 1/4" thick



### 6" ADJUSTABLE T-RAIL CLAMP

The adjustable T-rail clamp can be positioned along the grooved table top edge. Clamp is also removable.



#### 5" CASTER (OPTIONAL)

Recommended for mobile service providers.







717.759.3146 www.oakworksmed.com email: medsales@oakworks.com