



CFPM302 Imaging Table

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.



Integrated Headrest Table

Motions

- Motion 1 - Height Range 25" - 43"
- Motion 2 - Longitudinal Travel 10"
- Motion 3 - Lateral Tilt $\pm 15^\circ$

Features

- Choice of integrated headrest or rectangular carbon fiber tops with 56" of metal free imaging
- 2" Radiolucent Comfort Foam[®] table top pad
- Warranty is 3 years parts and 2 years labor

SPECIFICATIONS

- Large Adjustable Height Range 25" - 43"
- 56" Metal Free Imaging Space
- Cost Effective 3 Movement Table
- 500 lb Patient Weight Capacity
- Warranty: 3 years parts, 2 years labor

CFPM302 #84530 Integrated #84532 Rectangular

| INCLUDED FEATURES | |
|--|---|
| Carbon Fiber Top 56" Metal free imaging space | Integrated Headrest 22" (56cm) wide x 84" (214cm) long Rectangular 24" (61cm) wide x 80" (203cm) long |
| Motion 1 | Electric Height Range 25" - 43" (64cm-109cm) (without pad) |
| Motion 2 | Electric Longitudinal Travel 10" (25cm) |
| Motion 3 | Electric Lateral Tilt ±15° |
| Controls | All functions hand control |
| Weight Capacity | 500 lb (227kg) |
| Product Weight | 475 lb (216kg) |
| Included Features | 2" (5cm) Radiolucent Comfort Foam™ table top pad, 4" (10cm) locking casters, Crescent face rest pad (included with integrated headrest top), patient restraint straps (2) |
| Warranty | 3 years parts, 2 years labor |
| Safety Listings | FDA listed, CE marked, cETLus listed, CB Certificate |

OPTIONS must be ordered with table.

| | |
|----------------------------|--|
| 5" (13cm) casters | Recommended for mobile service providers |
| Foot Control | All functions control |
| Non-USA electrical systems | Please inquire |

ACCESSORIES order any time.

| ACCESSORIES | | # |
|-------------------------------------|-------------------------------------|-------|
| Arm Hammock | Use with Integrated Headrest Top | 48751 |
| Carbon Fiber Arm Board (6" wide) | 180° Articulation | 83575 |
| Fluoro Extender | Polycarbonate | 33975 |
| Interventional Supine Pillow System | Contoured pad and support wedge pad | 81821 |
| Spine Positioning System II | 8 piece system | 81819 |
| T-rail Clamp | 6" (15cm) long rail | 63336 |



Rectangular Top Table

ACCESSORIES

SPINE POSITIONING SYSTEM II

Radiolucency & adjustability for ideal cervical, thoracic & lumbar visualization. Minimize unwanted movement & reduce procedure time.



ARM HAMMOCK

This unique product will comfortably & easily support your patient's arms when prone. Easily attached to the Integrated Head Rest table top.



SINGLE CARBON FIBER ARM BOARD

Fully radiolucent with a 180° range of motion, the arm board requires no mounting or attachment hardware and simply slides under the table top pad. It can be used on any fluoroscopy table & will support arms and feet in a variety of positioning configurations.



FLUORO EXTENDER

Made of radiolucent, polycarbonate material, it's organic shape is conducive to a variety of table placement applications. 30" long x 23.5" wide x 1/4" thick



6" ADJUSTABLE T-RAIL CLAMP

The adjustable T-rail clamp can be positioned along the grooved table top edge. Clamp is also removable.



5" CASTER (OPTIONAL)

Recommended for mobile service providers.



717.759.3146 www.oakworksmid.com email: medsales@oakworksmid.com