# PT400M Physical Therapy Table

- Powered center break & height
- Optional arm rests & table extenders
- Contoured ergonomic design
- 550 lb. Weight Capacity









### PT400M

#PKG6099-T\* 1.75" Firm Response #PKG6770-T\* 2.5" Comfort Foam

FEATURES				
2 Section Contoured Top	29.5" (75cm) wide (at widest point) x 74" (188cm) long (including head rest)			
Motion I	Electric Height Range 18"-36" (44-90cm)			
Motion 2	Electric Center Break 0°-45°			
Controls	Multi-function foot control			
Weight Capacity	550 lb (250kg)			
Product Weight	186 lb (84kg)			
Padding	1.75"(4cm) Firm Response or 2.5"(6cm) Comfort Foam™			
Fabric	Premium, PVC free TerraTouch™			
Included Features	QuickLock™ Head Rest & Boiance™ pad, Accessory Mounting Plate (both ends)			
Warranty	3 years parts, 2 years labor			
Safety Listings	FDA listed, CE marked, cETLus approved, CB Certificate			

OPTIONS must be ordered with table		#
Locking Twin Casters (one end of the table)	To easily move your table in the room	63400
Non-USA electric system	call for ordering details	73145
Second Foot Control	Controls all powered movements	61900
Hand Control	Controls all powered movements	69405
ACCESSORIES order any time.		
Extension Pad	Use as table extender or lowered arm rest shelf (bracket required) 22.5" (57cm) wide x 9" (23cm) long	72187-T*
Rotating Arm Rests	set of 2 4.5"(11cm) wide x 16"(41cm) long	72184-T*
Shelf Bracket	to use Extension Pad as a lowered arm rest shelf	72203
Stabilization Strap	stabilize your patient	63530

#### STIFF TWIN BEAM CONSTRUCTION




Stiffer than cantilever designs and with fewer moving parts, our design supports all your work with a 18"- 36" height range and a powerful ETL listed motor.


#### FIRM RESPONSE FOAM

Firm Response foam is designed to deliver a firm yet comfortable surface that maximizes Physical Therapists' ergonomic efficiency during



examination and intervention procedures. Additionally, this surface minimizes Manipulation Dry Hits with its super high density, twin layer design. The Firm Response foam enables increased reliability of patient examinations.







**Dr. Robert DuVall,** PT, DHSc, ATC, OCS, SCS, FAAOMPT CSCS



## TERRATOUCH™ FABRIC

T26	T09	T21	T12	T18	T28	T07
Pure White	Butter	Saffron	Opal	Clay	Taupe	Earth
T27	T08	T19	T05	T20	T17	T04
Espresso	Ruby	Orchid	Heron	Ocean	Sapphire	Sky Blue
T30	T22	T24	T14	TI3	T29	T01
Seafoam	Sage	Blue Grass	Forest	Stone	Pewter	Coal

We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.

#### **CONTACT:**





Scan this code with your smart phone #72703

