

POWERLINE

Exam Table

- Cost effective solution for any facility.
- Partially assembled for quick set up.
- Hardwood construction boasts a 500 lb. weight capacity.
- 3 year parts warranty.



*Stone fabric, Grey finish
and storage shelf*



*Heron fabric, Natural Finish
and storage shelf*

TERRATOUGH™ FABRIC



We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.

FINISHES



PowerLine tables are highly water resistant and can be wiped down with disinfectants as needed.

Visit www.oakworksmed.com for a quote!



FM584901



WAVE BACKREST



The wave design makes adjusting your backrest quick and easy.

DURABLE MATERIALS



Hardwood construction, reinforced hardware, durable padding & fabric.

PAPER ROLL HOLDER



Paper Roll Holder is included standard.
Paper roll not included.

CONFIGURE YOUR TABLE

Backrest or Flat Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Height	30" (76cm)
Backrest Top	Manual Multi-Position
Weight Capacity	500 lb (227kg)
Product Weight	87 lb (39kg)
Padding	1.75" (4cm) Firm Response or 2.5" (6cm) Comfort Foam™
**Finish	Choice of: Natural and Grey
Paper Roll Holder	fits up to 22" (56cm) wide paper rolls
Storage Shelf	Fits both widths
Warranty	3 years parts, labor not included
Safety Listings	FDA listed, CE marked

STABILIZATION STRAP

Adjustable in length over a wide range for stabilization. Great for range of motion, mobilization, isolation.

**Watch
the video!**



Scan this code with your smart phone or visit the YouTube channel "Oakworks Medical"



ACCESSORIES <i>order any time.</i>		#
PowerLine Table Extender	27" (69cm) wide x 6" (15cm) long	51909-T*
PowerLine Table Extender	30" (76cm) wide x 6" (15cm) long	51910-T*
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support	PKG6914-T*
Stabilization Strap	Stabilize your patient	63530
Storage Shelf	Fits both widths	81576GR/81575NT
Prone Pillow	For prone or supine positioning	66180-T01

*Fabric color code must be specified.

** Finish color must be specified.

WEDGE BOLSTER

This product is not FDA listed.

Easily used as a backrest or leg support.
24" W x 10" H x 24" L with a 36° incline.



PRONE PILLOW

- Increased ventilation & improved comfort
- Low profile reduces cervical extension/flexion
- Superior positioning for any size patient
- Removable, easy to clean pads



717.759.3146 www.oakworksmed.com email: medsales@oakworks.com