



Task Light

LE-35

A better alternative to the multi-purpose gooseneck.



Philips Burton's new LE-35 provides two to five times more bright white light for your money than similar lights in its class. LE-35's K-arm reach of 42" (107 cm) is about double that of other comparable lights with drift-free and sag-free performance year after year. Our value, performance and two year product warranty make Philips Burton's LE-35 the best task light on the market.

The Right Light

PHILIPS
burton

Features

- Drift-free K-arm with 42" (107 cm) arm range
- Available in a table mount, wall mount and floorstand model
- UL 153 certified



LE-35 range of motion and general dimensions

Summary of technical data	
Illuminance:	30,000 lux (2787 fc) at 18" (46 cm)
Color Temperature:	2900 K
CRI:	97
Bulb life:	4000-hour average lamp life
Type of bulb:	Halogen Dichroic Reflector
Power:	35 Watts
Product warranty:	2-year warranty

Note: The optical data are nominal values based on measurements calculated according to UL/IEC performance standards at 120V. This product is NOT available in Norway, Sweden, Finland, Denmark, Estonia, Russia, Germany, UK, Ireland, Australia, Spain, the Netherlands or Iceland.

Ordering Information

120V Model	230/240V Model	Description
LE35FL	LE35FL25	Floorstand
LE35T	LE35T25	Table Mount
LE35W	LE35W25	Wall Mount
Accessories		
0009600PK		LE-35 Halogen Replacement Bulbs, 4/Box



Philips Burton
 21100 Lassen Street Chatsworth, CA 91311
 Phone: 818-701-8700 Fax: 818-701-8725
 www.burtonmedical.com
 © 2011 Philips Burton. All rights reserved. Product designs protected by copyright.
 We reserve the right to change details of design, materials and finishes; check website for most current product specifications.

LIT-0053, Rev. 3