



Magnifiers

Wave[®] Plus

A larger viewing area and 30% more intensity than conventional magnifier-type lights!


The Wave[®] Plus has a new sleek design to complement your office décor. Its unique, powerful 3.5 diopter magnifier is superior in both form and function. Plus, two independent fluorescent light bulbs offer 30% more intensity and enhance optical performance for your exams and surgical procedures.

The Right Light

PHILIPS
burton


Features

- Superior optics with extremely bright light source
- 41" (105 cm) Articulating K-Arm
- 6.5" x 4" (17 cm x 10 cm) lens
- "Quick-start" switch and ballast to activate bulbs
- New 'no drift' arm allows complete freedom of movement
- Available in floorstand and table/wall models



Summary of technical data	
Color temperature:	4100 K
Rated life of fluorescent bulb:	10,000 hours
Diopter of magnifying lens:	3.5
Focal length:	13" (33 cm)
Type of bulb:	2 independent 13-watt fluorescent bulb
Power requirements:	120V model standard; 230/240V models available
Certifications/Approvals:	UL 153 certified; CE marked
Product warranty:	1-year limited warranty

Note: The optical data are nominal values based on measurements calculated according to UL/IEC performance standards at 120V



Accessories

- | | |
|----------------|-------------------------------------------------|
| #0009700PK | Wave® Plus Replacement Bulbs, 120V, 4/Box |
| #0009310PK | Wave® Plus Replacement Bulbs, 230 / 240V, 4/Box |
| #06-800-140-00 | 4 Diopter Lens Attachment |
| #06-800-130-00 | 8 Diopter Lens Attachment |
| #06-800-120-00 | 10 Diopter Lens Attachment |



Philips Burton
 21100 Lassen Street Chatsworth, CA 91311
 Phone: 818-701-8700 Fax: 818-701-8725
 www.burtonmedical.com
 © 2011 Philips Burton. All rights reserved. Product designs protected by copyright.
 We reserve the right to change details of design, materials and finishes.

LIT-0008, Rev. 3