


SemperSure[®] Black

NITRILE EXAMINATION GLOVES



100 / 10 /

Skin-Friendly Black Examination Nitrile Gloves

- **Thickness:** Competes with a 3-4 mil glove
- Black color hides stains caused from blood, grease, dirt, dyes, and ink
- Manufactured without accelerators
- Low dermatitis potential
- Reduced risk of skin irritations
- Tested for use against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- ESD tested and approved
- Safe for use with sensitive electronics
- Meets FDA 21 CFR Part 177 requirements for food handling

REORDER INFO

Glove Color: Black

SKU	Size	Gloves/Box	Boxes/Case
SUNB101	X-Small	100 gloves/box	10 boxes/case
SUNB102	Small	100 gloves/box	10 boxes/case
SUNB103	Medium	100 gloves/box	10 boxes/case
SUNB104	Large	100 gloves/box	10 boxes/case
SUNB105	X-Large	100 gloves/box	10 boxes/case

SAMPLES: Available by request - ask your sales representative

For maximum protection, we recommend changing nitrile gloves every 20-30 minutes

For more information about protection against fentanyl, please contact your sales rep.

RECOMMENDED USES



PRODUCT SPECS



Rev07082022

sempermed[®]
YOUR GLOBAL GLOVE EXPERT™

100+ YEARS
- EST. 1920 -

Our Gloves are Disposable. The People they Protect are Not.®

www.sempermedusa.com • Ph 1.800.366.9545

