


sempermed®

supreme



quick don™

Latex Sterile Surgical Gloves

- **Thickness:** Competes with a 4-5 mil glove
- Surgical, 0.65 AQL
- Extra supple with great tactile sensitivity
- Fully anatomical with curved fingers and rolled cuff
- Excellent comfort and fit
- Micro-rough textured surface - superior instrument handling and grip
- Fast and easy to don: Quick Don™ inner coating
- Washed to minimize skin irritants
- Patented synthetic inner coating: US patent 6440499

RECOMMENDED USES



PRODUCT SPECS



REORDER INFO		Glove Color: Cream	
SPFP550	Size 5.5	50 pairs/box	6 boxes/case
SPFP600	Size 6	50 pairs/box	6 boxes/case
SPFP650	Size 6.5	50 pairs/box	6 boxes/case
SPFP700	Size 7	50 pairs/box	6 boxes/case
SPFP750	Size 7.5	50 pairs/box	6 boxes/case
SPFP800	Size 8	50 pairs/box	6 boxes/case
SPFP850	Size 8.5	50 pairs/box	6 boxes/case
SPFP900	Size 9	50 pairs/box	6 boxes/case

SAMPLES: Available by request - ask your sales representative

Rev12292021

