

sempermed® syntegra IR



quick don™

Polyisoprene Sterile Surgical Gloves

- **Thickness:** Competes with a 4-5 mil glove
- Surgical, 0.65 AQL
- Latex-free
- Fully anatomical with curved fingers and rolled cuff
- Similar molecular structure to natural rubber latex
- Manufactured with skin-friendly accelerators
- Low Dermatitis potential
- Fast and easy to don: Quick Don™ inner coating
- Micro-rough textured surface - superior instrument handling and grip

RECOMMENDED USES



PRODUCT SPECS



REORDER INFO		Glove Color: Cream	
SIR550	Size 5.5	50 pairs/box	6 boxes/case
SIR600	Size 6	50 pairs/box	6 boxes/case
SIR650	Size 6.5	50 pairs/box	6 boxes/case
SIR700	Size 7	50 pairs/box	6 boxes/case
SIR750	Size 7.5	50 pairs/box	6 boxes/case
SIR800	Size 8	50 pairs/box	6 boxes/case
SIR850	Size 8.5	50 pairs/box	6 boxes/case
SIR900	Size 9	50 pairs/box	6 boxes/case

SAMPLES: Available by request - ask your sales representative

Rev02142022

