

## User Manual

Item # PS058 - SleepSense® Female SNAP Adapter

**Thank you for purchasing a SleepSense® sensor.**

### Indication for Use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body Position or limb movements, for recording on an approved data acquisition system they are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

### Instructions

The cable is intended to connect two bands in series, to use on XXXL patients.

Any two 900X bands can be connected.

Use as in the pictures



## Warnings and Precautions

- SleepSense® sensors are for professional use only.
- Caution: US Federal Law restricts this device to sale, distribution and use by or on the order of a physician.
- SleepSense® sensors may be used only in conjunction with an approved recording system.
- Sensors are intended for use only on healthy skin.
- Discontinue use if product shows signs of wear, damage, or exposed metal

### **Limited Warranty**

SLP warrants the Inductive Interface cable to be free of defects in materials and workmanship for a period of twelve months from the date purchased. The sole liability of SLP and our Distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our Distributor(s) be liable for any loss of revenues or damage, direct, Consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This Warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to total. abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.



Manufactured by:

**S.L.P. Ltd.**  
62 Anilevitz St.  
Tel Aviv 67060  
ISRAEL  
Tel.: +972 3 6870640  
Fax: +972 3 5371282  
info@slp-med.com



S.L.P Inc.:

**S.L.P. Inc.**  
1121 E. Main St. #  
240, St. Charles, IL  
60174, USA  
Tel: +1 630 513 7479,  
Fax: +1 630 513 8478

EC	REP
----	-----

**Medes Limited**  
5 Beaumont Gate,  
Shenley Hill, Radlett,  
Herts WD7 7AR, England  
Tel/Fax: +44 19238 59810