



RLM Series LED



RLM - 14W Magnifier LED

Sunnex RLM-LED is a magnifying light with superior light intensity for detailed inspection. The magnification is 5 diopter (2.25X) and the light source is a 14W LED ring light with minimal heat emission.

The rugged and stable arm has adjustable joints with levers for easy and safe positioning. Can be mounted directly or with included C-clamp; mobile base is available as an added cost. As an option there is an add-on lens (2.25X) available. RLM-LED is pure in style, easy to keep clean and robust.

Applications

- Wound care
- Diagnostics
- Exam

Benefits

- Magnification 2.25 times
- Adjustable joints
- Adjustable lamp head
- Minimal heat emission
- Dimmable

Technical Specifications

Supply voltage	100-240V AC 50/60Hz, 1,8m cordset
Light source	14W LED ring light
Life span	> 25.000h
LED cover lens	Opal white
Total arm Length	880mm
Color	Gray
Color temperature	5500K
Light output	1250Lm / 10000Lx at 15cm
IP rating	IP 20
Protection Class	I
Base	Direct mount, c-clamp or mobile base

VERSION	LAMP	PART NO.
Direct Mount/C-Clamp	14W	FG RLM-LED-120V
Mobile	14W	ASB 32805

*Optional Add-on lens 2.25x

ACC RLM5XMAG