

# Extra Small Reusable Cup & Earclip Electrodes



**Technomed understands the importance of a crystal clear signal, ease of use and electrode durability.**

**Technomed's Reusable EEG Cup & Earclip Electrodes are designed to provide reliable signals and longtime ease of use.**

**These electrodes are designed for easy adjustment to patient contours to avoid electrode lift. The connector is optimized for easy headbox insertion and extraction.**


## **User Benefits**

- Extra small reusable cup & earclip electrodes
- Clear and reliable signal
- Short cup-to-cable connection to reduce lift risk
- Optimal patient comfort with low profile disk, soft and flexible anti-tangle lead wires
- Optimized connector design for easy headbox insertion and extraction
- Long electrode life due to watertight design and Kevlar reinforced wires

## Reusable Cup Electrodes

Materials	Cup Size	Cup Height	Wire Length	Part Number
Silver Plated (Ag)	6 mm	1.7 mm	1 m	TE/C14-434
Silver Plated (Ag)	6 mm	1.7 mm	1.5 m	TE/C14-634
Silver Plated (Ag)	6 mm	1.7 mm	2 m	TE/C14-834
Gold Plated (Au)	6 mm	1.7 mm	1 m	TE/C34-434
Gold Plated (Au)	6 mm	1.7 mm	1.5 m	TE/C34-634
Gold Plated (Au)	6 mm	1.7 mm	2 m	TE/C34-834
Silver/Silver Chloride (Ag/AgCl)	6 mm	2.3 mm	1.5 m	TE/C24-634/99


**Packing:** sets of 12, 2 x 6 Colors (red, green, white, yellow, blue, brown)



## Reusable Earclip Electrodes

Materials	Cup Size	Cup Height	Wire Length	Part Number
Silver Plated (Ag)	6 mm	1.7 mm	1 m	TE/E14-439
Silver Plated (Ag)	6 mm	1.7 mm	1.5 m	TE/E14-639
Gold Plated (Au)	6 mm	1.7 mm	1 m	TE/E34-439
Gold Plated (Au)	6 mm	1.7 mm	1.5 m	TE/E34-639
Silver/Silver Chloride (Ag/AgCl)	6 mm	1.7 mm	1 m	TE/E24-439
Silver/Silver Chloride (Ag/AgCl)	6 mm	1.7 mm	1.5 m	TE/E24-639

**Packing:** 1 pair




## Specifications

Lead-free	✓
Latex-free	✓
BPA-free	✓
DEHP < 0.1% by weight	✓
REACH & RoHS compliant	✓
Lead wire design	Tinsel wire
Lead wire core	Kevlar and Silver Plated Copper
Lead wire insulation	TPU
Spring loaded connector	Gold Plated Brass / PVC
Connector	DIN 42802

**Cleaning and disinfection:** see IFU for details

Your Technomed Business Partner:



USA
Rx ONLY



Taking you further. Step by step.

### Technomed

Amerikalaan 71  
6199 AE | Maastricht-Airport | The Netherlands  
T +31 43 408 68 68 | info@technomed.nl

### Technomed USA

Suite 4990 Constellation Drive  
White Bear Township | MN 55127 | United States  
T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl