


# Next Generation PowerSpin™ HXV Centrifuge

UNICO® is proud to present the PowerSpin™ HXV Centrifuge with 6 and 8 place horizontal or 12 place fixed angle rotor, variable speed up to 3500 RPM to insure a clean and complete separation. Quiet and dependable, with a maintenance free motor, it can be used to process blood, urine, semen and other materials. A single speed option also available.


The all digital front panel makes it easy to select speed and run time. Large LCD display shows status, time, speed (RPM) and RCF (g). Precise speed control enables you to comply with CLIA recommendations for specific specimens.



- Zero RPM Safety Latching System prevents lid from being opened while rotor is spinning
- Quiet running, only 50 decibels at 3400 RPM
- 99 minute digital timer, 1 sec increment

- 500-3500 RPM, 100 RPM increment
- Remembers last setting for quick repeat runs
- Stainless steel rotor chamber for easy cleaning
- Suction cup rubber feet eliminate vibration

UNICO Model	Item	Speed	Max. RCF	Rotor	Tube Size	Timer
PowerSpin™ HXV	C8302	3500 RPM	1777	8 place horizontal	2-15 ml	99 min digital
PowerSpin™ HXV	C8306	500 – 3500 RPM	1777	6 place horizontal	2-15 ml	99 min digital
PowerSpin™ HXV	C8308	500 – 3500 RPM	1777	8 place horizontal	2-15 ml	99 min digital
PowerSpin™ HXV	C8312	500 – 3500 RPM	1777	12 place fixed angle	2-15 ml	99 min digital



## UNICO Professional Digital Non-Contact Photo-Tachometer

- Highly accurate, scientific tool, easy to operate
- Aim and measure--RPMs instantly display on large LCD screen.
- Includes custom protective cushioned storage case and tape
- Order pn# C800-75

UNITED PRODUCTS & INSTRUMENTS

182 Ridge Road, Suite E Dayton, New Jersey 08810 Tel: 800-588-9776 \* Fax: 732-274-1155

E-mail: [info@unicosci.com](mailto:info@unicosci.com) \* Web: [www.unicosci.com](http://www.unicosci.com)

© 2015 UNICO®

# UNICO®