

# Independence At Home Mobile Centrifuge™

Designed for home based primary care teams to allow healthcare providers to spend more time with their patients and deliver care tailored to the individual needs and preferences.

Safe, reliable and easy to operate it can be run from the standard car outlet or inside with the optional AC adapter.

Built in level and steel carrying handle adds mobility and ease of operation.

Optional Tube Holdster™ - side mounted balance tube rack



C826H with Tube Holdster™

## Optional accessories for Porta-Spin C826

Part #	Description	UOM
C826-06	DC power adapter	each
C826-12	AC Adapter	each
C800-18	Balance Tube Holdster™	each
C800-03	Short Tube Adapter	Pack of 10
C800-04	Pediatric Tube Adapter	Pack of 2

UNICO Model	Item #	Tube Holdster™	Timer	Speed	Max. RCF (g)	Angled Rotor
Porta-Spin™	C826	No	Mechanical 30 min	3500 RPM	1640	6 place
	C828	No				8 place
	C826H	Yes				6 place
	C828H	Yes				8 place



United Products & Instruments, Inc  
[www.unicosci.com](http://www.unicosci.com) \* [info@unicosci.com](mailto:info@unicosci.com)

Independence At Home Mobile Centrifuge™, Tube Holdster™ are Registered Trademarks of United Products & Instruments, Inc.

182 Ridge Road, Suite E Dayton, NJ 08810  
 Tel: 732-274-1155 \* Fax: 732-274-1151