

POWER TABLES



Clinton takes power tables to new heights with a complete line of power exam, ultrasound and imaging tables designed and engineered to fit the needs of physicians, technicians and therapists.

Clinically proven in the field, Clinton power tables, are among the best in their class. With quality construction, dependable smooth, operation and enough power to satisfy nearly every need, Clinton Power Tables are practical and, best of all, affordable.

Welded steel construction makes them tough enough for years of worry-free service, and two inches of firm foam padding topped with premium, knit-backed vinyl makes patient comfort a priority in all our models. Choose from Shrouded Power Exam Tables, Open Base Power Exam Tables, Open Base Power Imaging Tables, Power Ultrasound and Imaging Trendelenburg Tables and our most economical power table yet, our E-Series, Power Tables.

Make Clinton your one source for affordable, convenient power tables — and plug into the power of Clinton.

INDEX

Shrouded, Power Tables2-5
 Open Base, Power Tables.....6-9
 Open Base, Power, Imaging Tables.....10-11
 Power, Ultrasound, Trendelenburg Tables....12-15
 Power, Imaging, Trendelenburg Tables.....16-19
 Power Table Options20
 Color Selector..... 21



MADE IN THE USA



Clinton power tables are among the most practical and affordable power exam tables on the market. They offer smooth dependable power, stylish design and years of dependable service.

Clockwise from left to right: 81360—Shrouded Power Table with Adjustable Backrest, Drop Section & Stirrups; 81330—Shrouded Power Table with Adjustable backrest & Drop Section; 80079—Multi-Use, Power Imaging, Trendelenburg Table with Adjustable Backrest, Drop Section & Stirrups; 80330 Open Base Power Imaging Table with Adjustable Backrest & Drop Section.

SHROUDED, POWER TABLES

Shrouded, Power Tables feature:

- Ultra-strong, uni-frame constructed base and top frame provides a rigid patient platform
- All welded steel base and top frame
- 3/4" plywood top for extra support (1.9 cm)
- Textured ABS plastic encloses lift mechanism and makes cleaning easier
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder is internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge



Shrouded, Power Table with One Piece Top

- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 81100 | length | width | height range |
|---------|-----------|----------|---------------------|
| overall | 72" | 27" | 21"-34" |
| | 182.87 cm | 68.58 cm | 53.34 cm – 86.36 cm |

Shrouded, Power Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 81200 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 27" | 21"-34" |
| | 182.87 cm | 68.58 cm | 53.34 cm – 86.36 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |



81200



81100



Shown with optional, solid drop section

81330

Shrouded, Power Table with Adjustable Backrest & Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (specify when ordering)
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 81330 | length | width | height range |
|--------------|----------|----------|---------------------|
| overall | 76" | 27" | 21"-34" |
| | 193 cm | 68.58 cm | 53.34 cm – 86.36 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| drop section | 13" | 27" | |
| | 33 cm | 68.58 cm | |

Shrouded, Space Saver Power Table with Drop Section

- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (*page 20*)

| 81220 | length | width | height range |
|--------------|--------|----------|---------------------|
| overall | 76" | 27" | 21"-34" |
| | 193 cm | 68.58 cm | 53.34 cm – 86.36 cm |
| drop section | 13" | 27" | |
| | 33 cm | 68.58 cm | |



81220

Shrouded, Power Adjustable, Casting Table with ClintonClean™ Leg Rest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- ClintonClean™ Solid Plastic Surface leg rest is adjustable from either side and raises up to 30° from flat position
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (*page 20*)

| 81350 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 21"-34" |
| | 193 cm | 68.58 cm | 53.34 cm – 86.36 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| leg rest | 13" | 26" | |
| | 33 cm | 66 cm | |



81350

Shrouded, Power Table with Stirrups, Adjustable Backrest & Footrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- Retractable, 2-position, adjustable stirrups
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (*page 20*)

| 81360 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 21"-34" |
| | 193 cm | 68.58 cm | 53.34 cm – 86.36 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 13" | 19" | |
| | 33 cm | 48.26 cm | |



81360

Shrouded Extra-Wide Bariatric Straight Power Table

- All welded, uni-frame steel base and top frame
- 550 lbs. load capacity under normal use (249.47 kg)
- Optional side rails, locking casters & more (*page 20*)

| 84108-34 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 34" | 21"-34" |
| | 182.87 cm | 86.36 cm | 53.34 cm – 86.36 cm |

84108-40 40" (101.60 cm) width model



84108-40

Shrouded, Extra-Wide, Bariatric, Power Table with Adjustable Backrest

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat
- 550 lbs. load capacity under normal use (249.47 kg)
- Optional side rails, locking casters & more (*page 20*)

| 84208-34 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 34" | 21"-34" |
| | 182.87 cm | 86.36 cm | 53.34 cm – 86.36 cm |
| backrest | 38" | 34" | |
| | 96.52 cm | 86.36 cm | |

84208-40 40" (101.60 cm) width model



84208-34

Shrouded, Extra-Wide, Bariatric, Power Table with Adjustable Backrest & Drop Section

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 550 lbs. load capacity under normal use (249.47 kg)
- Optional side rails, locking casters & more (*page 20*)

| 84438-34 | length | width | height range |
|--------------|----------|----------|---------------------|
| overall | 76" | 34" | 21"-34" |
| | 193 cm | 86.36 cm | 53.34 cm – 86.36 cm |
| backrest | 38" | 34" | |
| | 96.52 cm | 86.36 cm | |
| drop section | 13" | 34" | |
| | 33 cm | 86.36 cm | |

84438-40 40" (101.60 cm) width model



84438-34

OPEN BASE, POWER TABLES

Open Base, Power Tables feature:

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- 3/4" plywood top for extra support (1.9 cm)
- All welded steel base, top frame and scissors lift
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge

Open Base, Power Table with Stirrups, Adjustable Backrest & Footrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- Retractable, 2-position, adjustable stirrups
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 80360 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 13" | 19" | |
| | 33 cm | 48.26 cm | |



Open Base, Power Table with One Piece Top

- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 80100 | length | width | height range |
|-------|-----------|----------|---------------------|
| | 72" | 27" | 18"-35" |
| | 182.87 cm | 68.58 cm | 45.71 cm - 88.90 cm |



Open Base, Power Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 80200 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| Backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |



Open Base, Power Table with Adjustable Backrest & Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (specify when ordering)
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 80330 | length | width | height range |
|--------------|----------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| drop section | 13" | 27" | |
| | 33 cm | 68.58 cm | |



Space Saver Power Table with Drop Section

- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (*page 20*)

| 80220 | length | width | height range |
|--------------|--------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| drop section | 13" | 27" | |
| | 33 cm | 68.58 cm | |



80220

Power Adjustable Casting Table with ClintonClean™ Surface Leg Rest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- ClintonClean™ Solid Plastic Surface leg rest is adjustable from either side and raises up to 30° from flat position
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (*page 20*)

| 80350 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| leg rest | 13" | 26" | |
| | 33 cm | 66 cm | |



80350

Open Base, Extra-Wide, Bariatric, Straight Power Table

- All welded, uni-frame steel base and top frame
- 550 lbs. load capacity under normal use (249.47 kg)
- Optional side rails, locking casters & more (*page 20*)

| 84100-34 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 34" | 18"-35" |
| | 182.87 cm | 86.36 cm | 45.71 cm - 88.90 cm |

84100-40 Same as above with 40" width (101.60 cm)



84100-40

Open Base, Power Table with Adjustable Backrest

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat
- All welded, uni-frame steel base and top frame
- 550 lbs. load capacity under normal use (249.47 kg)
- Optional side rails, locking casters & more (*page 20*)

| 84200-34 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 34" | 18"-35" |
| | 182.87 cm | 86.36 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 34" | |
| | 96.52 cm | 86.36 cm | |

84200-40 Same as above with 40" widths (101.60 cm)



84200-40

Open Base, Power Table with Adjustable Backrest & Drop Section

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 550 lbs. load capacity under normal use (249.47 kg)
- Optional side rails, locking casters & more (*page 20*)

| 84430-34 | length | width | height range |
|--------------|----------|----------|---------------------|
| overall | 76" | 34" | 18"-35" |
| | 193 cm | 86.36 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 34" | |
| | 96.52 cm | 86.36 cm | |
| drop section | 13" | 34" | |
| | 33 cm | 86.36 cm | |

84430-40 Same as above with 40" widths (101.60 cm)



84430-40

OPEN BASE, POWER IMAGING TABLES

Open Base, Power Imaging Tables feature:

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- 3/4" plywood top for extra support (1.9 cm)
- All welded steel base, top frame and scissors lift
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism that does not need to roll on the floor
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge

Open Base, Multi-Use, Power Imaging Table with Stirrups

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- Retractable, 2-position, adjustable stirrups
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 85309 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 13" | 19" | |
| | 33 cm | 48.26 cm | |



85309

Flat position with window drop engaged

Open Base, Power Imaging Table with Window Drop

- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 85100 | length | width | height range |
|-------|-----------|----------|---------------------|
| | 72" | 27" | 18"-35" |
| | 182.87 cm | 68.58 cm | 45.71 cm - 88.90 cm |



85100

Open Base, Power Imaging Table with Adjustable Backrest and Window Drop

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base, top frame and scissors lift
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 85200 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 27" | 18"-35" |
| | 182.87 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |



85200

Power Imaging Table with 3 Section Top

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- 450 lbs. load capacity under normal use (204 kg)
- Optional side rails, locking casters & more (page 20)

| 85300 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 18"-35" |
| | 193 cm | 68.58 cm | 45.71 cm - 88.90 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 13" | 27" | |
| | 33 cm | 68.58 cm | |



85300

**POWER, ULTRASOUND,
TREDELENBURG TABLES**

Power, Ultrasound, Trendelenburg Tables feature:

- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- 600 lbs load capacity under normal use (272 kg)



General, Flat Top, Ultrasound Table

- One piece flat top with firm foam padding
- Optional side rails, locking casters & more (page 20)

| 80061 | length | width | height range |
|-------|-----------|----------|---|
| | 72" | 27" | 21 ¹ / ₂ " – 37 ¹ / ₂ " |
| | 182.88 cm | 68.58 cm | 54.61 cm – 95.25 cm |



80061



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

General Ultrasound Table with Three-Section Top

- Three-piece top with pneumatic adjustable backrest and manual footrest
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails, locking casters & more (page 20)

| 80063 | length | width | height range |
|----------|----------|----------|---|
| overall | 76" | 27" | 21 ¹ / ₂ " – 37 ¹ / ₂ " |
| | 193 cm | 68.58 cm | 54.61 cm – 95.25 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 13" | 27" | |
| | 33 cm | 68.58 cm | |



80063



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

Multi-Use, Ultrasound Table with Stirrups

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails, locking casters & more (page 20)

| 80069 | length | width | height range |
|----------|----------|----------|---|
| overall | 76" | 27" | 21 ¹ / ₂ " – 37 ¹ / ₂ " |
| | 193 cm | 68.58 cm | 54.61 cm – 95.25 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 12" | 18.5" | |
| | 30.48 cm | 47 cm | |



80069



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

POWER, ULTRASOUND,
TRENDLENBURG TABLES

**General, Ultrasound Table with
Adjustable Backrest**

- Two-piece top with easy-to-adjust, pneumatic backrest
- Optional side rails, locking casters & more (page 20)

80062

| | | | |
|---------|-----------|----------|---|
| overall | 72" | 27" | 21 ¹ / ₂ " – 37 ¹ / ₂ " |
| | 182.88 cm | 68.58 cm | 54.61 cm – 95.25 cm |

| | | | |
|----------|----------|----------|------------------------------------|
| backrest | 38" | 27" | POWER Backrest Available |
| | 96.52 cm | 68.58 cm | |



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters



Power Imaging, Trendelenburg Tables feature:

- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- 600 lbs load capacity under normal use (272 kg)



Flat Top, Power Imaging Table with Drop Window

- One piece flat top with drop window
- Optional side rails, locking casters & more (page 20)

| 80071 | length | width | height range |
|---------|-----------|----------|---------------------|
| overall | 72" | 27" | 21 1/2" – 37 1/2" |
| | 182.88 cm | 68.58 cm | 54.61 cm – 95.25 cm |
| window | 10" | 10" | |
| | 25.4 cm | 25.4 cm | |



80061



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

Power Imaging Table with Fowler Back and Drop Window

- Two-piece top with pneumatic adjustable backrest and window drop
- Optional side rails, locking casters & more (page 20)

| 80072 | length | width | height range |
|----------|-----------|----------|---------------------|
| overall | 72" | 27" | 21 1/2" – 37 1/2" |
| | 182.88 cm | 68.58 cm | 54.61 cm – 95.25 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| window | 10" | 10" | |
| | 25.4 cm | 25.4 cm | |



80072



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

Power Imaging Table with Three-Section Top and Drop Window

- Three-piece top with pneumatic adjustable backrest, manual footrest and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails, locking casters & more (page 20)

| 80073 | length | width | height range |
|----------|----------|----------|---------------------|
| overall | 76" | 27" | 21 1/2" – 37 1/2" |
| | 193 cm | 68.58 cm | 54.61 cm – 95.25 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 12" | 27" | |
| | 30.48 cm | 68.58 cm | |
| window | 10" | 10" | |
| | 25.4 cm | 25.4 cm | |



80073



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

POWER, IMAGING,
TRENDLENBURG TABLES

Multi-Use, Power Imaging Table with Stirrups and Drop Window

- Three-section table with pneumatic adjustable backrest, manual footrest, stirrups and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails, locking casters & more (page 20)

| 80079 | length | width | height range |
|----------|----------|----------------------------------|---|
| overall | 76" | 27" | 21 ¹ / ₂ " – 37 ¹ / ₂ " |
| | 193 cm | 68.58 cm | 54.61 cm – 95.25 cm |
| backrest | 38" | 27" | |
| | 96.52 cm | 68.58 cm | |
| footrest | 12" | 18 ¹ / ₂ " | |
| | 30.28 cm | 47 cm | |
| window | 10" | 10" | |
| | 25.4 cm | 25.4 cm | |

POWER
Backrest & Footrest
Available



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters



5 Year
Limited
Warranty

POWER TABLE OPTIONS

Safety Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set

| 080 | length | width |
|-----|-----------|-------|
| | 92" | 2" |
| | 233.67 cm | 5 cm |

Safety Rails (pair)

- 3 height position
- Plunger height adjustment

082

2-Hook IV Pole and Mounting Bracket

- Chrome pole with height adjustment
- Removable bracket mounts to either side of head end
- Also available with 4-hook top (must specify)

083

Power Backrest

(Trendelenburg Tables only)

- Smooth power control raises backrest to 80° and lowers to flat

011

Power Footrest

(Trendelenburg Tables only)

- Smooth power control lowers footrest to 30° and raises to flat

090

Hand-Held Height Control

(Not available on E-Series Tables)

- Alternative to foot control
- Stores easily under top when not in use

084

Locking Casters

(Not available on E-Series Tables)

- Heavy duty, rubber wheel, easy-roll casters

087

30" Width Top

(Trendelenburg Tables only)

037



COLOR SELECTOR


UPHOLSTERY

Clinton uses tough, woven, knit-backed, vinyl upholstery that resists bunching and stretching. Its durable surface can be sanitized quickly and easily with your specific disinfectant and a soft cloth.

(See detailed cleaning instructions with product.)

In an effort to be current on changing color trends and to better provide our customers with the most up-to-date palette, please visit our web site at www.clinton-ind.com.

A selection of Clinton upholstery colors available for each specific product can be found on every product page by clicking the "Color Chart" link.





800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

www.clinton-ind.com