



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Three Drawer Infection Control Cart with Breakaway Lock, Standard Package #6500



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Stainless steel break-away lock bar
- Pack of 50 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Full wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Aluminum side rail / push handle
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled
- Twelve year warranty

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Isolation Cart Upgrade:

- Silver ion antimicrobial paint available

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Assembled Cart Dimensions

- 41.5" H x 32" W x 22" D
(105.4 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 115 pounds (52.2 kg)
- Approx. shipping weight: 140 pounds (72.7 kg)