



EasyOne Pro / Pro LAB

Complete point of care pulmonary function testing

Portable DLCO, lung volumes & spirometry.
A highly efficient and accurate way to test lung function.



- EasyOne Pro**
- Spirometry
 - DLCO
 - Provocation testing
- Lung volumes**
- Single breath helium dilution

EasyOne Pro LAB

- Spirometry
 - DLCO
 - Provocation testing
- Lung volumes**
- Multi-breath nitrogen washout



Simple. Certain. Proven.

Features

- Automated user guidance & test quality feedback based on 2019 & 2005 ATS/ERS standards
- Z-score, LLN and % predicted for fast interpretation of results
- Real-time curves and pediatric incentives
- Export of pdf files and raw data
- Flexible HL7 and XML interface for easy EMR integration
- Hygienic solution with Spirette and Barriette consumables eliminates the risk of cross-contamination
- Compact device with smooth surfaces for easy and thorough cleaning

Specifications

EasyOne Pro: FVC, FVL, SVC, MVV, Pre/Post, Provocation, DLCO

EasyOne Pro LAB: FVC, FVL SVC, MVV, Pre/Post, Provocation, DLCO, MBW, LCI

Languages: English, Chinese, Croatian, Danish, Dutch, Finnish, French, German, Italian, Japanese, Norwegian, Portuguese, Russian, Spanish, Swedish, Turkish, Vietnamese

EMR connectivity: HL7, XML, GDT, API with EasyOne Connect

Dimensions: 10.6 x 13.2 x 10.6" / 17.65lb / 27 x 33.5cm³ (HXWXD), 8kg

Part numbers

3000-1: EasyOne Pro

3100-1: EasyOne Pro LAB

2050-1: Spirette box of 50

2050-5: Spirette box of 200

2050-10: Spirette box of 500

3050-1: DLCO Barriettes box of 50

3050-2: DLCO Barriettes box of 100

3150-1: FRC Barriette box of 40 (for use with EasyOne Pro LAB only)

3150-2: FRC Barriette box of 80 (for use with EasyOne Pro LAB only)

TrueFlow ^{ndd} makes the difference

The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity.

ndd TrueFlow is a hygienic and resistance-free solution that does not require calibration during its lifetime.

EasyOne Connect

EasyOne Connect software offers one database and one platform for all EasyOne point-of-care-solutions. Data management has never been easier.

TrueCheck ^{ndd}

TrueCheck takes care of the essential quality control for gas analysis testing. EasyOne Pro is the only device proven to be accurate for a lifetime for DLCO measurements.

