

A photograph of a person's leg with a prosthetic knee joint submerged in clear blue water. The prosthetic is black and blue, with a silver metal shaft. The person is wearing a light blue patterned swimsuit and a white prosthetic foot. The water is clear with some bubbles and light reflections.

HIGH-FUNCTIONING & VERSATILE KNEES FOR ALL MOBILITY LEVELS

Exceptionally stable on all terrains.

Suitable for use in the shower, pool, at the beach,
and for other water activities.

Promotes confidence and proper gait mechanics
that result in a more natural, energy-efficient stride.

Now with threaded tops for lower knee clearance!

For more information visit legworks.com.

sales@legworks.com | (408) 692-5633

ALL-TERRAIN KNEE



ATK-PA-01 (PYRAMID ADAPTOR)
ATK-TT-01 (THREADED TOP)

FRESH WATERPROOF
ATK-SF-01 (STANCE FLEXION
 with PYRAMID ADAPTOR)

MOBILITY LEVEL	K1-K4		
WEIGHT RATING	275 lb (125 kg)		
FLEXION ANGLE	150°		
DISTAL CONNECTION	30 mm Tube Clamp		
AMPUTATION LEVEL	Transfemoral, Knee Disarticulation		Transfemoral
PRODUCT HEIGHT	PA: 7 5/16 in (186 mm)	TT: 7 1/16 in (179 mm)	7 15/16 in (201 mm)
PRODUCT WEIGHT	PA: 2.29 lb (1,040 g)	TT: 2.30 lb (1,045 g)	2.71 lb (1,230 g)
SUGGESTED CODES*	L5611, L5850, L5925		L5611, L5845, L5850, L5925

ALL-TERRAIN KNEE PREMIUM



ATK-PA-02 (PYRAMID ADAPTOR)
ATK-TT-02 (THREADED TOP)

WATERPROOF & CORROSION RESISTANT
ATK-SF-02 (STANCE FLEXION
 with PYRAMID ADAPTOR)

MOBILITY LEVEL	K1-K4		
WEIGHT RATING	330 lb (150 kg)		
FLEXION ANGLE	150°		
DISTAL CONNECTION	30 mm Tube Clamp		
AMPUTATION LEVEL	Transfemoral, Knee Disarticulation		Transfemoral
PRODUCT HEIGHT	PA: 7 5/16 in (186 mm)	TT: 7 1/16 in (179 mm)	7 15/16 in (201 mm)
PRODUCT WEIGHT	PA: 2.30 lb (1,044 g)	TT: 2.31 lb (1,049 g)	2.72 lb (1,234 g)
SUGGESTED CODES*	L5611, L5850, L5925		L5611, L5845, L5850, L5925

ATK HYDRAPRO



HYD-PA-03 (PYRAMID ADAPTOR)
HYD-TT-03 (THREADED TOP)

WATERPROOF & CORROSION RESISTANT HYDRAULIC
HYD-SF-03 (STANCE FLEXION
 with PYRAMID ADAPTOR)

MOBILITY LEVEL	K3-K4		
WEIGHT RATING	330 lb (150 kg)		
FLEXION ANGLE	150°		
DISTAL CONNECTION	30 mm Tube Clamp		
AMPUTATION LEVEL	Transfemoral, Knee Disarticulation		Transfemoral
PRODUCT HEIGHT	PA: 7 5/16 in (186 mm)	TT: 7 1/16 in (179 mm)	7 15/16 in (201 mm)
PRODUCT WEIGHT	PA: 2.32 lb (1,052 g)	TT: 2.33 lb (1,057 g)	2.74 lb (1,242 g)
SUGGESTED CODES*	L5814, L5850, L5925		L5814, L5845, L5850, L5925

*L-Codes suggested by LegWorks are for reference-only, and are subject to change. The responsibility and liability for coding lies with the provider.