



**LegWorks**

## **HD SERIES**

# **HIGH-FUNCTIONING & VERSATILE KNEES FOR ALL MOBILITY LEVELS**

***NOW WITH A 440 LB WEIGHT RATING***

Exceptionally stable on all terrains.

Suitable for use in the shower, pool, at the beach,  
and for other water activities.

Promotes confidence and proper gait mechanics  
that result in a more natural, energy-efficient stride.

**For more information** visit [legworks.com](http://legworks.com).

[sales@legworks.com](mailto:sales@legworks.com) | (408) 692-5633



## WATERPROOF & CORROSION RESISTANT

### ALL-TERRAIN KNEE PREMIUM HD

ATK-PA-03 (PYRAMID ADAPTOR)  
ATK-TT-03 (THREADED TOP)

ATK-SF-03 (STANCE FLEXION with  
PYRAMID ADAPTOR)



<b>MOBILITY LEVEL</b>	K1-K4		
<b>WEIGHT RATING</b>	440 lb (200 kg)		
<b>FLEXION ANGLE</b>	150°		
<b>DISTAL CONNECTION</b>	30 mm Tube Clamp		
<b>AMPUTATION LEVEL</b>	Transfemoral, Knee Disarticulation		Transfemoral
<b>PRODUCT HEIGHT</b>	PA: 7 5/16 in (186 mm)	TT: 7 1/16 in (179 mm)	7 15/16 in (201 mm)
<b>PRODUCT WEIGHT</b>	PA: 2.30 lb (1,044 g)	TT: 2.31 lb (1,049 g)	2.72 lb (1,234 g)
<b>SUGGESTED CODES*</b>	L5611, L5850, L5925		L5611, L5845, L5850, L5925



## WATERPROOF & CORROSION RESISTANT HYDRAULIC

### ATK HYDRAPRO HD

HYD-PA-03 (PYRAMID ADAPTOR)  
HYD-TT-03 (THREADED TOP)

HYD-SF-03 (STANCE FLEXION with  
PYRAMID ADAPTOR)



<b>MOBILITY LEVEL</b>	K3-K4		
<b>WEIGHT RATING</b>	440 lb (200 kg)		
<b>FLEXION ANGLE</b>	150°		
<b>DISTAL CONNECTION</b>	30 mm Tube Clamp		
<b>AMPUTATION LEVEL</b>	Transfemoral, Knee Disarticulation		Transfemoral
<b>PRODUCT HEIGHT</b>	PA: 7 5/16 in (186 mm)	TT: 7 1/16 in (179 mm)	7 15/16 in (201 mm)
<b>PRODUCT WEIGHT</b>	PA: 2.32 lb (1,052 g)	TT: 2.33 lb (1,057 g)	2.74 lb (1,242 g)
<b>SUGGESTED CODES*</b>	L5814, L5850, L5925		L5814, L5845, L5850, L5925

\*L-Codes suggested by LegWorks are for reference-only, and are subject to change. The responsibility and liability for coding lies with the provider.