



ATK
HYDRAPRO
with optional stance flexion

Your everyday waterproof hydraulic knee

FEATURES & BENEFITS



HYDRAULIC SWING

The patent-pending Hydraulic Variable Cadence Controller smoothly facilitates multiple walking speeds on all terrain.



CORROSION RESISTANT

A waterproof prosthesis, the ATK HydraPro is suitable for use in the shower, pool, at the beach, and for other water activities.



RUGGED & DURABLE

The ATK HydraPro's advanced composite design is weight rated to 330 lb (150 kg) and comes with a 24-month warranty.



STABLE & EFFICIENT

Patented AutoLock Technology promotes confidence and proper gait mechanics that result in a more natural, energy-efficient stride.



INDICATED USERS

Active amputees seeking a hydraulic, waterproof primary prosthesis that is low-maintenance.

Amputees seeking a secondary prosthesis with a waterproof and corrosion resistant hydraulic knee (e.g. shower/pool/beach leg).

Amputees that desire a highly stable hydraulic knee.

ASK ABOUT OUR 30-DAY LOANER KNEE PROGRAM!



ATK HYDRAPRO HYD-PA-02

ATK HYDRAPRO WITH STANCE FLEXION HYD-SF-02

MOBILITY LEVEL	K3-K4	
WEIGHT RATING	330 lb (150 kg)	
FLEXION ANGLE	150°	
PROXIMAL CONNECTION	Pyramid Adaptor	
DISTAL CONNECTION	30 mm Tube Clamp	
AMPUTATION LEVEL	Transfemoral, Knee Disarticulation	Transfemoral
PRODUCT HEIGHT	7 5/16 in (186 mm)	7 15/16 in (201 mm)
PRODUCT WEIGHT	2.32 lb (1,052 g)	2.74 lb (1,242 g)
SUGGESTED CODES	L5814, L5850, L5925	L5814, L5845, L5850, L5925

*L-Codes suggested by LegWorks are for reference-only, and are subject to change. The responsibility and liability for coding lies with the provider.