

ORIGINAL SERIES

# ALL-TERRAIN KNEE

with optional Stance Flexion



## SIMPLIFIED STABILITY

A simple twist-knob engages the manual lock mode, providing additional safety in early rehabilitation or situational use.



## PURELY MECHANICAL

The Variable Cadence Controller facilitates multiple walking speeds without the use of hydraulics or pneumatics.



## STABLE & EFFICIENT

AutoLock Technology promotes proper gait mechanics that result in a more natural, energy efficient stride.



## 2-YEAR GUARANTEE

24-month limited manufacturer's warranty

# ALL-TERRAIN KNEE

ORIGINAL SERIES

## FEATURES & BENEFITS

- Stability provided via the patented AutoLock Technology which engages at full swing phase extension
- Encourages K2 amputees to transition to the K3 mobility level via the patent-pending Variable Cadence Controller
- Increased safety on uneven terrain when simple twist-knob manual lock is engaged
- Durable and weight rated to 275 lb (125 kg)
- Waterproof in fresh water (eg. shower, lake)



## INDICATED USERS

- K2 amputees with potential to progress to K3 mobility
- Underinsured amputees who would benefit from a high-functioning and cost-effective prosthetic knee
- Amputees who would benefit from a knee with inherent stability, including hip disarticulation and bilateral
- Amputees in preparatory phase of prosthetic delivery

## SUGGESTED L-CODES\*

L5611, L5845<sup>†</sup>, L5850, L5925

<sup>†</sup>Available on stance flexion models only

## WARRANTY

24-Month limited manufacturer's warranty

## ALL-TERRAIN KNEE

### ATK-PA-01



<b>AMPUTATION LEVEL</b>	Transfemoral, Knee Disarticulation	<b>FLEXION ANGLE</b>	150°
<b>MOBILITY LEVEL</b>	K1-K4	<b>PROXIMAL CONNECTION</b>	Pyramid Adaptor
<b>WEIGHT RATING</b>	275 lb (125 kg)	<b>DISTAL CONNECTION</b>	30 mm Tube Clamp
<b>PRODUCT WEIGHT</b>	2.29 lb (1,040 g)	<b>CLEARANCE</b>	1 <sup>5</sup> / <sub>16</sub> in (24 mm)
<b>PRODUCT HEIGHT</b>	7 <sup>5</sup> / <sub>16</sub> in (186 mm)	<b>FITTED HEIGHT</b>	4 <sup>15</sup> / <sub>16</sub> in (125 mm)

## ALL-TERRAIN KNEE WITH STANCE FLEXION

### ATK-SF-01



<b>AMPUTATION LEVEL</b>	Transfemoral	<b>FLEXION ANGLE</b>	150°
<b>MOBILITY LEVEL</b>	K1-K4	<b>PROXIMAL CONNECTION</b>	Pyramid Adaptor
<b>WEIGHT RATING</b>	275 lb (125 kg)	<b>DISTAL CONNECTION</b>	30 mm Tube Clamp
<b>PRODUCT WEIGHT</b>	2.71 lb (1,230 g)	<b>CLEARANCE</b>	1 <sup>9</sup> / <sub>16</sub> in (40 mm)
<b>PRODUCT HEIGHT</b>	7 <sup>15</sup> / <sub>16</sub> in (201 mm)	<b>FITTED HEIGHT</b>	5 <sup>5</sup> / <sub>16</sub> in (141 mm)

\*L-Codes suggested by LegWorks are for reference-only, and are subject to change. The responsibility and liability for coding lies with the provider.