

**PREMIUM SERIES**

# ALL-TERRAIN KNEE

with optional Stance Flexion



## CORROSION RESISTANT

An ideal secondary waterproof prosthesis, the All-Terrain Knee Premium is great for the pool, beach, and other water activities.



## RUGGED & DURABLE

The All-Terrain Knee Premium's advanced composite design is weight rated to 330 lb (150 kg) and comes with a 24-month limited manufacturer's warranty.



## PURELY MECHANICAL

The Variable Cadence Controller facilitates multiple walking speeds without the use of hydraulics or pneumatics.



## STABLE & EFFICIENT

AutoLock Technology promotes proper gait mechanics that result in a more natural, energy efficient stride.

# ALL-TERRAIN KNEE

PREMIUM SERIES

## FEATURES & BENEFITS

- Stability provided via the patented AutoLock Technology which engages at full swing phase extension
- Worry-free variable cadence without the use of hydraulics or pneumatics
- Waterproof and corrosion resistant (for use in the shower, pool, ocean, and other outdoor activities)
- Durable and weight rated to 330 lb (150 kg)
- Situational use of manual lock allows for stability in or around water



## INDICATED USERS

- Amputees seeking a high-functioning, low-maintenance knee for everyday use
- Amputees seeking a secondary waterproof prosthesis
- K2 amputees with potential to progress to K3 mobility
- Underinsured amputees who would benefit from a high-functioning and cost-effective prosthetic knee

## SUGGESTED L-CODES\*

L5611, L5845<sup>†</sup>, L5850, L5925

<sup>†</sup>Available on stance flexion models only

## WARRANTY

24-Month limited manufacturer's warranty

## ALL-TERRAIN KNEE PREMIUM

ATK-PA-02



<b>AMPUTATION LEVEL</b>	Transfemoral, Knee Disarticulation	<b>FLEXION ANGLE</b>	150°
<b>MOBILITY LEVEL</b>	K1-K4	<b>PROXIMAL CONNECTION</b>	Pyramid Adaptor
<b>WEIGHT RATING</b>	330 lb (150 kg)	<b>DISTAL CONNECTION</b>	30 mm Tube Clamp
<b>PRODUCT WEIGHT</b>	2.30 lb (1,044 g)	<b>CLEARANCE</b>	1 <sup>5</sup> / <sub>16</sub> in (24 mm)
<b>PRODUCT HEIGHT</b>	7 <sup>5</sup> / <sub>16</sub> in (186 mm)	<b>FITTED HEIGHT</b>	4 <sup>1</sup> / <sub>16</sub> in (125 mm)

## ALL-TERRAIN KNEE PREMIUM WITH STANCE FLEXION

ATK-SF-02



<b>AMPUTATION LEVEL</b>	Transfemoral	<b>FLEXION ANGLE</b>	150°
<b>MOBILITY LEVEL</b>	K1-K4	<b>PROXIMAL CONNECTION</b>	Pyramid Adaptor
<b>WEIGHT RATING</b>	330 lb (150 kg)	<b>DISTAL CONNECTION</b>	30 mm Tube Clamp
<b>PRODUCT WEIGHT</b>	2.72 lb (1,234 g)	<b>CLEARANCE</b>	1 <sup>9</sup> / <sub>16</sub> in (40 mm)
<b>PRODUCT HEIGHT</b>	7 <sup>15</sup> / <sub>16</sub> in (201 mm)	<b>FITTED HEIGHT</b>	5 <sup>5</sup> / <sub>16</sub> in (141 mm)

\*L-Codes suggested by LegWorks are for reference-only, and are subject to change. The responsibility and liability for coding lies with the provider.