

EZwaste[®] Heavy Duty Solvent Waste System

www.foxxlifesciences.com



MADE

ISO 9001 • ISO 13485

life cycle may vary



Our Innovation. Your Opportunity.

EZwaste[®] Heavy Duty Solvent Waste System

www.foxxlifesciences.com




CONTAINER SIZE	1/8" OD PORTS	1/4" OD PORTS	1/4" or 3/8" HB PORT	PART NUMBER
2.5 Liter	6			332-1412-OEM
4 Liter	6			332-9412-OEM
5 Liter	6			332-2412-OEM
10 Liter	6			332-3412-OEM
20 Liter	6			332-4412-OEM
20 Liter Round	6			332-D412-OEM
2.5 Liter	6		1	332-1418-OEM
4 Liter	6		1	332-9418-OEM
5 Liter	6		1	332-2418-OEM
10 Liter	6		1	332-3418-OEM
20 Liter	6		1	332-4418-OEM
20 Liter Round	6		1	332-D418-OEM
2.5 Liter	4	3		332-1413-OEM
4 Liter	4	3		332-9413-OEM
5 Liter	4	3		332-2413-OEM
10 Liter	4	3		332-3413-OEM
20 Liter	4	3		332-4413-OEM
20 Liter Round	4	3		332-D413-OEM
2.5 Liter	4	4		332-1421-OEM
4 Liter	4	4		332-9421-OEM
5 Liter	4	4		332-2421-OEM
10 Liter	4	4		332-3421-OEM
20 Liter	4	4		332-4421-OEM
20 Liter Round	4	4		332-D421-OEM
2.5 Liter	4	3	1	332-1419-OEM
4 Liter	4	3	1	332-9419-OEM
5 Liter	4	3	1	332-2419-OEM
10 Liter	4	3	1	332-3419-OEM
20 Liter	4	3	1	332-4419-OEM
20 Liter Round	4	3	1	332-D419-OEM


* 1/16" OD Adapters are also avaliable for the HD line as cap only assemblies.

Solvent Waste System Sizing Guide

Due to all the possibilities with configuring your EZwaste Solvent Waste System, we have provided the following images to help identify the fitting required for your laboratory setup. All measurements are approximate.



2 For measuring ID, place end of tube flat on the piece of paper and determine if the dark area is contained by the tubing.



MADE

Rev. 33-001-02 01/28/2016

Foxx Life Sciences, 6 Delaware Drive, Salem, NH 03079 www.foxxlifesciences.com Order Now: 1 (603) 890-FOXX (3699)