

THE WORLD'S **ONLY** DOUBLE FDA CLEARED MICRONEEDLING DEVICE

FOR PROVEN WRINKLE & ACNE SCAR REDUCTION



DISCOVER OUR **NEW**MICRONEEDLING DEVICE

EXCEED

AMIEA MED EXCEED

The medical microneedling device with proven results.

With our effective NEW device you can perform percutaneous collagen induction (PCI) on your patient's skin. A minimally invasive proceedure, microneedling uses a needle to puncture the epidermis and dermis. The process ensures the skin does not suffer open wounds and epithelial function is restored quickly and the treatment requires only a very short regeneration phase.

CLINICAL PERFORMANCE

To property assess the clinical performance of our microneedling system a clinical evaluation trial has been conducted. The study has proven that four microneedling treatments, in intervals of 4 weeks, lead to measurable skin improvement. Furthermore a smoother and rejuvenated skin appearance was achieved. The treatment was very well tolerated with minimal downtime, giving effective client satisfaction.

HIGHLIGHTS & KEY FACTS

PRECISE PERFORMANCE THROUGH POWERFUL HANDPIECE

- · No tearing of tissue through strong penetration technician and retraction with optimized drive
- · Regulated penetration depth through tilting needle head with 90° angle penetration of the skin
- · Small cartridge head offers possibility to easily treat areas around the eyes, nose and mouth

HIGHEST STANDARDS FOR QUALITY AND SAFETY

- · Scale for needle length indication to allow precise depth control
- · Safety membrane in the disposable cartridge prevents contamination of the handpiece
- · Specific safety features; e.g., device recognises if handpiece is not attached properly
- · Detachable cable at the handpiece to minimize breakage and failure

COMFORTABLE OPERATION

- · Optimal ergonomics for working without fatigue
- · Digital control unit
- · Compact Design
- · Selectable footswitch mode between continuous operation and ON/OFF





ACCESSORIES

- · Safety needle cartridges (sterile single-use products)
 - · 6-needleplate 1.5 Exceed
 - · 1-Liner Exceed
 - · 3-Liner Exceed
- foot pedal
- · Handpiece cover

CERTIFIED AREA OF APPLICATION

The treatment of acne scars on the face for the clinical improvement of scar appearance and in particular the reduction of scar depth and normalization of the skin structure on the chosen area of application.


CERTIFICATIONS AND PATENTS

- · Medical device certification (CE0123) according to MDD and DIN EN ISO 13485 standard
- · US 6,345,553; US 6,505,530; US 8,029,527
- · EP 1 495 782; EP 1 618 915; EP 1 872 823; EP 1 958 659; EP 2 462 979
- Additional patents pending
- FDA approval pending
- · Engineered and manufactured in Germany

FEATURES

BENEFITS

| Adjustable needle speed | High power results in deeper penetrationLess pain and minimal downtime |
|---------------------------------|---|
| Needle length freely selectable | Scale of 0 to 1,9 ensures reproducability of needle length setting Appropriate needle length for each area with our exchangeable cartridges Sensitive or aggressive treatment |
| Precision and power | Low epidermal traumaLow noise and vibrationSuperior results |
| No ablation of the epidermis | Minimal downtime Treatment of Fitzpatrick I–V Thickening of the epidermis |
| Specials | High-contrast digital displaySetting of foot-switch functionKey tone settings |





Unique Skin Ltd info@uniqueskin.co.uk +44 (0)1462 481500