



The Premier Freeze Drying System

The ultimate in freeze drying technology

If your lab is looking for the most modern, hassle-free freeze drying technology available for ice water hash, look no further than our Labconco FDry-8L freeze dryer system. The Labconco FDry-8L offers the ultimate in freeze drying technology with real-time monitoring and the ability to easily customize runs with programmable shelf temperatures and vacuum levels.

3x daily runs compared to other freeze dryers

8 hour production cycles allow your lab to complete three full runs each day, which is less than half the time of other freeze dryers.

Smart failsafes can save you \$30,000+ if something goes wrong

The Labconco's ethernet connectivity and built-in sensors allow for remote 24/7 remote monitoring and alerts. In the event of a problem, your team can react quickly to prevent any batch losses, which means you have peace of mind for every run.

Large throughput capacity

With 5x 10.25" x 19.50" trays, totaling 7.2 sq ft of shelf space, the Labconco FDry-8L freeze dryer is capable of drying up to 700 dry grams per run producing 2,100 dry grams per 24 hours.

Adjustable vacuum pressure

Unlike other freeze dryers, operators can fully calibrate their vacuum pressure to ensure supremely consistent batch drying and recipe fine tuning.

Smart controls

A vibrant 5" color touch screen with an intuitive user interface gives you intelligent and precise control of your freeze dryer's operation.

ETL & CE certified

Comes with industry standard electrical certifications that regulators and inspectors expect.



Labconco FDry-8L **\$18,995** + shipping

Includes premium pharmaceutical grade **Leybold SOGEVAC NEO 16 D vacuum pump** with purchase price.



1 year manufacturer's warranty

Available in 115v or 230v configurations, and requires a dedicated 20amp circuit for operation.



| [gopurepressure](https://www.gopurepressure.com)

[purepressure.com](https://www.purepressure.com)



Leybold SOGEVAC NEO D 16 vacuum pump

Vacuum pump maintenance as infrequently as every 3 years.

Internal baffling system prevents any oil aerosolization or back flow into the freeze dryer chamber, a common problem with cheaper oil pumps.

Capable of pulling vacuum down to .1 millibar.

Ultra quiet operation, only 54 decibels.

Contact our solventless equipment experts today!

sales@gopurepressure.com | 720-446-9565



GET PERFECT BATCH SETTINGS RIGHT OUT OF THE BOX WITH OUR EXPERT MADE CUSTOM FDRY-8L SPECIFIC SOP

AVAILABLE FOR \$750

LEARN MORE AT [PURECANNALABS.COM](https://purecannalabs.com)