




FACE SHIELD DRAWINGS


Original plans and assembly instructions at:
<https://www.protohaven.org/proto-shield/>

Plans have been modified to accept 1/2" to 1" wide elastic and use 8 1/2 x 11 transparency film as shield

Headband components made from .060 PETG



Shield made from .020 PETG or overhead transparency film



ELASTIC
1/2" to 1" Wide
x 15.5 inches Long



Proto Shield by Mat Thorne and Devin Montgomery is licensed under a Creative Commons Attribution-ShareAlike 4.0 International License. Based on a work at <https://www.prusaprinters.org/prints/25857>.

LEGAL: This design and these instructions are provided for free and as-is. Neither Protohaven nor the creators of the design and instructions offer any warranty or assume any liability for their use.