



Drug Detection And Identification Products

PRODUCTS AND CAPABILITIES

Drug Aerosols test technology detects both trace and substance level of drugs. Testing is convenient and efficient, available in multiple sizes.

Set of single and multi-drug ampoule based field test kits for the detection and identification of illicit drugs.

PDT is a series of individual presumptive reagent drug tests designed to quickly and easily identify controlled substances while in the field.



Drug Type	Aerosols	PenTests	Presumptive Drug Tests (PDT)
Amphetamines, Cannabis, Narcotics, Opiates	D4D	D4D	N/A
Amphetamines, Narcotic Alkaloids	D4D	D4D	Mayer's; Marquis
Amphetamines	D4D	D4D	Mandeline; Fentanyl
Barbiturates	Barbitusol	Barbitusol	Dille-Koppanyi
Bath Salts	D4D	D4D	MDVP; Mephedrone
Buprenorphine	D4D	D4D	N/A
Cocaine / Crack	Coca	Coca	Scott's
Ephedrine	D4D	D4D; Ephedrine	Chen's
Fentanyl	Coming Soon	Coming Soon	Fentanyl
GHB / GBL	N/A	N/A	GHB
Heroin: Black Tar	N/A	N/A	Mecke's
Heroin: White, Brown	D4D; Herosol	D4D; Herosol	Mecke's
Ketamine	D4D	D4D	Morris
LSD	D4D; Herosol	D4D; LSD	Ehrlich's
Marijuana	D4D; Cannabis	D4D; Cannabis	Duquenois-Levine; KN
Meth / Ecstasy	D4D; Meth-X	D4D; Meth-X	Methamphetamine; Ecstasy; Fentanyl
Methadone	D4D; Herosol	D4D	Methadone
Methaqualone	N/A	N/A	Methaqualone; PCP
Oxycodone; Codeine; Heroin, Morphine	D4D; Herosol	D4D; Opiatest	Opiates
PCP	D4D; Coca	D4D	PCP
Pentazocine	D4D	D4D	Talwin
PMA	D4D	D4D	Lieberman's
Propoxyphene	D4D	D4D	Propoxyphene
Psilocybin	D4D	D4D	Psilocybin
Spice / K2	D4D	D4D	N/A
Vallium / Rohypnol	N/A	BZO	Valium / Rohypnol