



**HI-LO TABLE MOTOR REPLACEMENT INSTRUCTIONS  
PN: 03105, MOTOR MAX65 8000N 120V, VENDOR MODEL  
MAX65-A200415A0100-000**

**CAUTION: MAKE CERTAIN MOTOR IS UN-PLUGGED AND TOP FRAME IS IN THE COMPLETE DOWN POSITION OR IS LYING ON ITS' SIDE BEFORE STARTING REPAIRS.**

Remove power and control switch cords as required using the special tool enclosed.

To replace motor, remove bolts, nuts and washers from each end of defective motor. If table top is elevated it should be lowered at this time to allow the top frame to rest on the base frame. Make certain power and control switch cords are securely attached to replacement motor before installing. The motor shaft may be turned in or out to lengthen or shorten the shaft if needed.