



## 美利馳 A1 電動輔助式後輪



## Merits A1 Power Pack for Manual Wheelchair

### Merits Power Pack for Manual Wheelchair



#### Specifications:

<b>max.speed</b>	5.5 kph / 3.4 mph
<b>range up</b>	16 km / 10 mi.
<b>motor</b>	DC 12V / 150W
<b>battery</b>	12V / 20AH
<b>charger</b>	1.5 A off-board
<b>drive wheel</b>	8" PU tire x 2 pcs
<b>capacity</b>	136 kg / 300 lb
<b>dimension</b>	300 L x 200 W x 245 H
<b>weight-base</b>	6.4 kg / 14 lb
<b>weight-battery</b>	6.8 kg / 15 lb
<b>weight-total</b>	13.2 kg / 29 lb
<b>wheelchair width</b>	14 in ~ 20 in

**Twin-wheel design with 136kg / 300lb weight capacity and much greater traction on loose gravel and wet grass**

**Easy to attach to any wheelchair, and simply walk behind with your thumb on the control and let the power pack do all the work for you**

**With speeds up to 5.5kph / 3.4mph and a range of up to 16km /10miles, there is more than enough capacity for most days out**

**Light enough to take anywhere, packs neatly into a travel bag enabling that extra mobility anywhere**