









| Switch                  | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..  | Feature                                       |
|-------------------------|---|------------|--------------------|------------------|--|--------------|--|---|
| <b>Big Red</b>          |    | Mechanical | 12.7 cm            | 156 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | All purpose and classroom use too.                             | Large activation area, 4 interchangeable caps |
| <b>Jelly Bean</b>       |    | Mechanical | 6.4 cm             | 71 g             | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | All purpose, classroom use.                                    | 4 interchangeable caps                        |
| <b>Spec</b>             |   | Mechanical | 3.5 cm             | 100 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | All purpose but smaller target area, great for individual use. | Small and discreet activation area            |
| <b>Big Buddy Button</b> |  | Mechanical | 11.5 cm            | 150 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | All purpose and classroom too.                                 | Large activation area, flush with base.       |





| Switch                | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..                     | Feature                                     |
|-----------------------|---|------------|--------------------|------------------|--|--------------|---------------------------------|---|
| <b>Buddy Button</b>   |    | Mechanical | 6.3 cm             | 100 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | All purpose and classroom too.  | Activation area flush with base.            |
| <b>Microlight</b>     |    | Mechanical | 4.5 x 1.3 cm       | 10 g             | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | Discreet personal use.          | Small, finger size.                         |
| <b>Big Candy Corn</b> |   | Proximity  | 9 cm               | 0 g              | <input checked="" type="checkbox"/> Auditory<br><input checked="" type="checkbox"/> Visual<br><input type="checkbox"/> Tactile | Ablenet      | Small hand or finger movements. | Big, touch free activation area.            |
| <b>Candy Corn</b>     |  | Proximity  | 4.7 cm             | 0 g              | <input checked="" type="checkbox"/> Auditory<br><input checked="" type="checkbox"/> Visual<br><input type="checkbox"/> Tactile | Ablenet      | Small hand or finger movements. | Small discreet, touch free activation area. |

| Switch          | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..   | Feature  |
|-----------------|---|------------|--------------------|------------------|--|--------------|---|--|
| <b>Pillow</b>   |    | Mechanical | 8 cm               | 180 g            | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile            | Ablenet      | Face activation.  | Soft, washable, activation area.                         |
| <b>Mini Cup</b> |    | Mechanical | 2.5 cm             | 130 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | Discreet personal use.  | Small discreet, activation area                          |
| <b>Grasp</b>    |   | Mechanical | 7.6 x 3.5 cm       | 300 g            | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile            | Ablenet      |   | Squeeze to activate.                                     |
| <b>Plate</b>    |  | Mechanical | 8 x 8 cm           | 114 g            | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile                       | Ablenet      | Making a book interactive and other creative accessibility hacks maximising the discreet form factor. | Discreet, flat form factor, activate with a light touch. |

| Switch                      | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..   | Feature                                |
|-----------------------------|---|------------|--------------------|------------------|--|--------------|---|--|
| <b>Wobble</b>               |    | Mechanical | 10 cm              | 132 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | Users who have a waving, striking motion, may reduce discomfort or possible pain when impacting with a surface without that much give, like a button. | Strike from any direction to activate. |
| <b>Blue 2</b>               |    | Mechanical | 7 x7 cm            |                  | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | iPads, Chromebooks and PC computers.  | Wireless, Bluetooth                    |
| <b>Piko Small WR</b>        |   | Mechanical | 3 cm               | 125 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Piko         | Discreet size, for personal use outdoors like on a wheel chair.   | Rugged, water proof                    |
| <b>Piko Water Resistant</b> |  | Mechanical | 5 cm               | 125 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Piko         | Outdoors, kitchen or personal use on wheelchair.  | rugged, water proof                    |



| Switch               | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..   | Feature  |
|----------------------|---|------------|--------------------|------------------|--|--------------|---|--|
| <b>Piko</b>          |    | Mechanical | 3 cm/5 cm          | 75 g / 125 g     | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Piko         | Heavy use and high traffic areas in a classroom or individual use                                     | Ruggedized.  |
| <b>iSwitch</b>       |    | Mechanical | 7.5 cm             | 120 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Inclusive    | Classroom use wirelessly. Pass the switch around to run whole class projected activities.             | Wireless with large activation area.                     |
| <b>Pal Pad Mini</b>  |   | Mechanical | 5 x 5 cm           | 57 g             | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile                       | Adaptivation | Making a book interactive and other creative accessibility hacks maximising the discreet form factor. | Discreet, flat form factor, activate with a light touch. |
| <b>Pal Pad Small</b> |  | Mechanical | 7 x 10 cm          | 57 g             | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile                       | Adaptivation | Making a book interactive and other creative accessibility hacks maximising the discreet form factor. | Discreet, flat form factor, activate with a light touch. |





| Switch                 | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer       | Great for..   | Feature  |
|------------------------|---|------------|--------------------|------------------|--|--------------------|---|--|
| <b>Pal Pad Large</b>   |    | Mechanical | 11.5 x 15.25 cm    | 57 g             | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile | Adaptivation       | Making a book interactive and other creative accessibility hacks maximising the discreet form factor. | Discreet, flat form factor, activate with a light touch. |
| <b>Smoothie</b>        |    | Mechanical | 12.5 cm            | 120 g            | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile | Inclusive          | All purpose and classroom use too.  | Large activation area, flush with base.                  |
| <b>Ultra Light 1HD</b> |   | Mechanical |                    |                  | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile | Marblesoft         | Classroom, or individual use for discreet and light touch.  | Finger size, stick on anywhere                           |
| <b>Orby</b>            |  | Mechanical | 6.35 cm            | 100 g            | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile | Origin Instruments |   | Wireless   |

| Switch              | Image   | Type       | Activation Surface | Activation force | Feedback  | Manufacturer | Great for..  | Feature   |
|---------------------|---|------------|--------------------|------------------|---|--------------|--|---|
| <b>Cosmo</b>        |    | Mechanical | 5 cm               | 50 to 500 g      | <input checked="" type="checkbox"/> Auditory<br><input checked="" type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Filisia      | Classroom, shared activities, and beginning/reluctant switch users. Great for turn taking and participation. | Wireless, Adjustable activation and audio and light feedback. |
| <b>Honeybee</b>     |    | Proximity  | 14 x 7 cm          | 0 g              | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile                       | Adaptivation | Small hand or finger movements.  | Adjustable sensing range.                                     |
| <b>Jelly Beamer</b> |   | Mechanical | 6.35 cm            | 71 g             | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile            | Ablenet      | Classroom or individual wireless, use. Pass the switch around to run whole class projected activities.       | Wireless  |
| <b>Easy Switch</b>  |  | Mechanical | 7.5 cm             |                  | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile            | Inclusive    | Easy setup of wireless switch activities on PC's with available USB port.                                    | Wireless  |

| Switch                    | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..  | Feature  |
|---------------------------|---|------------|--------------------|------------------|--|--------------|--|--|
| <b>Bright Box Tactile</b> |    | Mechanical | 10.7 x 6.3         |                  | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile                       | Adaptivation | Great for beginning/reluctant switch users who need extra feedback.                                    | Lights up for additional feedback.                                       |
| <b>Big Beamer</b>         |    | Mechanical | 12.7 cm            | 156 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Ablenet      | Classic profile with wireless access. Pass around the switch for shared, whole class activities.       | Wireless   |
| <b>Table Tapper</b>       |   | Mechanical |                    |                  | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile                       | Adaptivation | Creative accessibility hacks maximising the discreet form but large activation area.                   | Turns a table top, counter or other surface into a tap-sensitive switch. |
| <b>SimplyWorks Large</b>  |  | Mechanical | 12.5 cm            | 120 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Inclusive    | Classroom or individual wireless, use. Pass the switch around to run whole class projected activities. | Wireless, low profile  |



| Switch                      | Image   | Type       | Activation Surface | Activation force | Feedback   | Manufacturer | Great for..  | Feature                                       |
|-----------------------------|---|------------|--------------------|------------------|--|--------------|--|---|
| <b>SimplyWorks Small</b>    |    | Mechanical | 7.5 cm             | 120 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Inclusive    | Classroom or individual wireless, use, in a smaller footprint. Pass the switch around to run whole class projected activities.   | Wireless, low profile                         |
| <b>Pneumatic Switch Kit</b> |    | Mechanical |                    |                  | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile            | Adaptation   | New or reluctant switch users who can grasp or squeeze.  | Turns squeeze toys into switches.             |
| <b>Piko Dual Small</b>      |   | Mechanical | 3 cm               | 125 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Piko         | User who needs two switches but close enough for one finger activation. Rugged construction perfect for wheelchair use.  | Ruggedized, double switch in small footprint. |
| <b>Piko Dual Medium</b>     |  | Mechanical | 5 cm               | 125 g            | <input checked="" type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Piko         | User who needs two switches but close enough for one finger activation. Rugged construction perfect for wheelchair use. Larger activation for multi-users in classroom, too. | Ruggedized, double switch in small footprint. |
| <b>Piko FingerButton</b>    |  | Mechanical | 1 cm               | 75 g             | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile            | Piko         | User with minimal finger movements   | Light touch, velcros around finger.           |

| Switch                  | Image   | Type       | Activation Surface | Activation force | Feedback  | Manufacturer | Great for..  | Feature   |
|-------------------------|---|------------|--------------------|------------------|---|--------------|--|---|
| <b>Piko Finger Dual</b> |    | Mechanical | 1 cm               | 75 g             | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input checked="" type="checkbox"/> Tactile | Piko         | User with minimal finger movements                                 | light touch   |
| <b>Taction Pads</b>     |    | Conductive |                    |                  | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile            | Adaptavation | Flexible positioning, for user with minimal movement and dexterity | Attaches to virtually any surface, activated by skin contact. Connect with powered Taction Pad Adapter. |
| <b>Piko 5 function</b>  |   | Mechanical | 1.7 x 4.9 cm       | 75/125 g         | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile            | Piko         | Computer control with a single finger                              | 5 activation areas in a very small footprint. Requires Hitch interface                                  |
| <b>Flexible Switch</b>  |  | Mechanical | 10 cm, 7.5 x 13 cm | 57 g             | <input type="checkbox"/> Auditory<br><input type="checkbox"/> Visual<br><input type="checkbox"/> Tactile            | Adaptavation | Flexible positioning, for user with minimal movement and dexterity | Attaches to virtually any surface, activated by skin contact. Connect with powered Taction Pad Adapter. |

**Visual:** lights up      **Tactile:** feel the switch moving

**Auditory:** beep, click or tone